





**ΔΙΑΤΜΗΜΑΤΙΚΟ ΠΡΟΓΡΑΜΜΑ ΜΕΤΑΠΤΥΧΙΑΚΩΝ ΣΠΟΥΔΩΝ
«ΕΦΑΡΜΟΓΕΣ ΤΗΣ ΒΙΟΛΟΓΙΑΣ ΣΤΗΝ ΙΑΤΡΙΚΗ»
Μάθημα: Μικροβιολογία και Δημόσια Υγεία**



Μέθοδοι διάγνωσης της μικροβιακής αντοχής

Ολυμπία Ζαρκωτού,

Διευθύντρια Μικροβιολογικού Τμήματος, Γ.Ν. Πειραιά «Τζάνειο»






AMR
a "slow-moving
tsunami"

Antibiotic resistance: a growing threat to human health

Antibiotic resistance is the ability of bacteria to combat the action of one or more antibiotics. Bacteria, not humans or animals, become antibiotic-resistant.

In Europe, the health impact of antibiotic resistant infections is comparable to that of influenza, tuberculosis and HIV/AIDS combined.




More than
35000 deaths

Each year, more than 35 000 people die from antibiotic-resistant infections in the European Union, Iceland and Norway. This is equivalent to the number of passengers on 13 cruise ships.

Antibiotic resistance is a silent pandemic and a growing threat to human health.

13 cruise ships



Over **70 %**
healthcare-associated
infections

Over 70% of the health impact of antibiotic-resistant infections is directly linked to healthcare-associated infections. This could be minimized through adequate infection prevention and control measures, as well as antibiotic stewardship in healthcare settings.

Increasing burden


Resistance to antibiotics that are used as last line for treatment of infections, such as the carbapenems, has the highest health impact.

Between 2016 and 2020, the overall number of deaths caused by antibiotic-resistant bacteria under study has increased.

For carbapenem-resistant *Klebsiella pneumoniae* and *Acinetobacter* spp, commonly causing healthcare-associated infections, the number of attributable deaths increased by approximately 50%.



Global deaths (counts) attributable to bacterial AMR by pathogen–drug combination, 2019




Antimicrobial Resistance Collaborators. Global burden of bacterial antimicrobial resistance in 2019: a systematic analysis. Lancet. 2022;399(10325):629-655.

Carbapenem-resistant *Acinetobacter baumannii*


Balkan Peninsula



Percentage of isolates with resistance



D Carbapenem-resistant *Acinetobacter baumannii*
Raw data



Carbapenem-resistant *Klebsiella pneumoniae*

Balkan Peninsula




Percentage of isolates with resistance




F Carbapenem-resistant *Klebsiella pneumoniae*

Raw data




Global burden of bacterial antimicrobial resistance in 2019: a systematic analysis. Lancet, 2022, 399: 629-55

Estimations of the burden of infections with antibiotic-resistant bacteria presented as attributable deaths per 100 000 population by country*, EU/EEA, 2020



Assessing the health burden of infections with antibiotic-resistant bacteria in the EU/EEA, 2016-2020. Stockholm: ECDC; 2022.



Στόχοι μικροβιακής αντοχής: Ελλάδα



EUROPEAN
ANTIBIOTIC
AWARENESS DAY

A EUROPEAN
HEALTH INITIATIVE



ecdc

EUROPEAN CENTRE FOR
DISEASE PREVENTION AND
CONTROL


Στόχοι μικροβιακής αντοχής: πώς μπορούμε να τους επιτύχουμε έως το 2030;



<https://antibiotic.ecdc.europa.eu/en/publications-data/country-sheet-2024-greece>

Antimicrobial resistance targets ⁱ -2024 update ⁱⁱ - (2023 data)			
Greece			
	Target achieved	Progress	Regress
 Reduce by 27% the total consumption of antibiotics in humans Defined daily doses (DDDs) per 1 000 inhabitants per day	2019 baseline	34.1	-
	2023	28.5	-16.3%
	2030 TARGET	24.9	-27%
 At least 65% of the total consumption of antibiotics in humans belongs to the 'Access' group of antibiotics As defined in the AWaRe classification of the WHO *Percentage point difference from 2019.	2019 baseline	46.8%	-
	2023	42.0%	-4.8%*
	2030 TARGET	65%	+18.2%*
 Reduce by 10% the total incidence of bloodstream infections with meticillin-resistant <i>Staphylococcus aureus</i> (MRSA) Number per 100 000 population	2019 baseline	4.6	-
	2023	6.5	+41.8%
	2030 TARGET	4.1	-10%
 Maintain at baseline level the total incidence of bloodstream infections with third-generation cephalosporin-resistant <i>Escherichia coli</i> Number per 100 000 population	2019 baseline	2.6	-
	2023	5.6	+117%
	2030 TARGET	2.6	-
 Reduce by 5% the total incidence of bloodstream infections with carbapenem-resistant <i>Klebsiella pneumoniae</i> Number per 100 000 population	2019 baseline	13.1	-
	2023	21.4	+64.3%
	2030 TARGET	12.4	-5%

Εκτιμώμενη συνολική επίπτωση βακτηριαιμιών με ανθεκτικό φαινότυπο (η ανά 100 000 πληθυσμού), 2019-2023, Ελλάδα



Μεταβολή %, 2019-2023, ανά παθογόνο

- Στελέχη *Escherichia coli* με αντοχή στις κεφαλοσπορίνες τρίτης γενιάς: αύξηση 117,4%
- Ανθεκτικά στις καρβαπενέμες στελέχη *K. pneumoniae*: αύξηση 64,3%
- Ανθεκτικά στις καρβαπενέμες στελέχη *P. aeruginosa*: αύξηση 97,2%
- Ανθεκτικά στις καρβαπενέμες στελέχη *A. baumannii*: αύξηση 63,2%
- MRSA: αύξηση 41,8%
- *E. faecium* VRE: αύξηση 166%

European Centre for Disease Prevention and Control.
Antimicrobial resistance in the EU/EEA (EARS-Net)- Annual
Epidemiological Report 2023. Stockholm: ECDC; 2024.

TACKLING ANTIMICROBIAL RESISTANCE ON TEN FRONTS



Public awareness



Sanitation and hygiene



Antibiotics in agriculture and the environment



Vaccines and alternatives



Surveillance



Rapid diagnostics



Human capital



Drugs



Global Innovation Fund



International coalition for action

Review on Antimicrobial Resistance

'Tackling drug-resistant infections globally: Final report and recommendations'. The review on antimicrobial resistance; May 2016



A European One Health Action Plan against Antimicrobial Resistance (AMR)

https://health.ec.europa.eu/system/files/2020-01/amr_2017_action-plan_0.pdf

Σύσταση του Συμβουλίου σχετικά με την ενίσχυση των δράσεων της ΕΕ για την καταπολέμηση της μικροβιακής αντοχής στο πλαίσιο της προσέγγισης «Μία υγεία» 2023/C 220/01

[https://eur-lex.europa.eu/legal_content/EL/TXT/?uri=CELEX:32023H0622\(01\)](https://eur-lex.europa.eu/legal_content/EL/TXT/?uri=CELEX:32023H0622(01))



EN

Search

Research and innovation

Home > ... > All research and innovation news > New European Partnership on One Health AMR: €253 million for research and innovation against antimicrobial resistance

NEWS ARTICLE | 23 September 2025 | Directorate-General for Research and Innovation | 2 min read

New European Partnership on One Health AMR: €253 million for research and innovation against antimicrobial resistance

Today, the European Union has officially launched the [Partnership on One Health Antimicrobial Resistance \(OHAMR\)](#). Co-funded by [Horizon Europe](#) with €75 million, this ambitious 10-year programme tackles the urgent global health threat of antimicrobial resistance (AMR).

European Partnership on One Health Antimicrobial Resistance

Kick-off

23 September 2025
Brussels, Belgium

OHAMR

Co-funded by the European Union

The European Partnership on One Health AMR (EUP OHAMR)

brings together 53 partners from 30 countries to strengthen AMR R&I with a One Health approach.

By coordinating research, boosting innovation, improving surveillance and **diagnostics**, and supporting better prevention strategies, the partnership aims to close critical knowledge gaps and translate science into policy

<https://cordis.europa.eu/project/id/101217154>

‘Best buy’ interventions to tackle AMR



- improving hygiene in healthcare facilities, including promotion of handwashing
- stewardship programmes promoting more prudent use of antibiotics to end decades of over-prescription

These will decrease the health burden of AMR – measured in DALYs – by about 40%

- **use of rapid diagnostic tests**
- delayed prescription
- public awareness campaigns


Source: OECD (2018), Stemming the Superbug Tide: Just a Few Dollars More. Available at: oe.cd/amr-2018.


Μικροβιολογικό εργαστήριο και ανίχνευση αντοχής




Πληροφορία







β- λακταμάσες: υδρολυτικό φάσμα



*including benzylpenicillins, aminopenicillins, carboxypenicillins, ureidopenicillin, narrow spectrum cephalosporins (cefazolin and cefuroxime and others)

**BSBL substrates plus oxacillin, nafcillin, and dicloxacillin

β- λακταμάσες και νεότερα αντιβιοτικά


	ESBL	KPC	MBL	AmpC	OXA-48	<i>P. aeruginosa</i> (MDR/XDR)	<i>Acinetobacter</i> (MDR/XDR)	<i>S. maltophilia</i>
Aztreonam/avibactam								
Cefepime/enmetazobactam								
Cefepime/taniborbactam								
Cefepime/zidebactam								
Cefiderocol								
Ceftaroline/avibactam								
Ceftolozane/tazobactam								
Ceftazidime/avibactam								
Imipenem/relebactam								
Meropenem/nacubactam								
Meropenem/vaborbactam								

Μέθοδοι ανίχνευσης της αντοχής

Phenotypic methods

*Constrained by
bacterial growth time*

Dilution



48 h

Diffusion




48 h

Gradient test



48 h

Chromogenic media




24 h

Automated devices



> 20 h

LFIA




Molecular-based methods

PCR Gene 1 Gene 2... Gene N




several h

qPCR



Cycle sequencing



24 - 48 h

NGS




several days

** estimated time depend on a sample analysed
(clinical specimen vs. isolated bacterial culture)*

Mass spectrometry

MALDI-TOF MS

Incubation (antibiotic + sample)



Matrix



several h


** estimated time depend on a sample analysed
(clinical specimen vs. isolated bacterial culture)*

Table 1. Differences between phenotypic and genotypic resistance detection methods

Characteristic	Phenotypic methods	Genotypic methods
Question to be answered	Does the antibiotic inhibit bacterial growth at clinically relevant concentrations?	Is a gene or mutation associated with antibiotic resistance present?
Turn-around time	Slow	Fast 
Inoculum needed	High	Low
Provides information about resistance mechanism	No	Yes
Predicts antibiotic susceptibility and resistance	Yes	Sometimes. Only detects a gene or mutation associated with resistance; this may not correlate with phenotypic resistance in all isolates (e.g. if a gene is not expressed). If a singular genotypic resistance type is associated with resistance to a particular antibiotic in a particular bacterial species, its absence infers susceptibility. However, when there is more than one genotypic resistance type associated with resistance, such an inference may not always be correct.
Provides MIC	Yes	No
Cost	Moderate 	High




Κλινικό
δείγμα




Ζώνες
αναστολής

ID / AST



Τιμές MIC




Pheno R Μηχανισμός;

TAT= 48-72 h



CIDT



TAT= 1-4 h




ID / Ανίχνευση γονιδίων αντοχής

Διάγνωση
Κατάλληλη
αγωγή



AST



Ταχύτητα

αξιοπιστία

επαναληψιμότητα



Προαγωγή της φροντίδας
του ασθενούς



Διασφάλιση συνταγαγρόφησης
του κατάλληλου AB




Προαγωγή της ορθολογικής
χρήσης AB



Ανίχνευση μηχανισμών αντοχής
Εφαρμογή μέτρων για περιορισμό
της διασποράς

Φαινοτυπικές μέθοδοι: Μέθοδος διάχυσης των δίσκων (Kirby-Bauer)



Bauer A.W., Kirby W.M.M., Sherris J.C. and Turck M. 1966.

Antibiotic susceptibility testing by a standardized single disk method. Am. J. Clin. Pathol. 45:493-496

EUCAST disk diffusion method for AST. Version 13.0, January 2025.

Έλεγχος ευαισθησίας και ερμηνεία των αποτελεσμάτων



Μη μέλη: \$200.00

Επιδημιολογικά cut-offs
Ενδογενής αντοχή
Ανίχνευση μηχανισμών αντοχής
Δοσολογικά σχήματα για SDD όρια

HomeContactSitemapNewsletter

EUCAST

EUROPEAN COMMITTEE
ON ANTIMICROBIAL
SUSCEPTIBILITY TESTING

European Society of Clinical Microbiology and Infectious Diseases

QUICK NAVIGATION

EUCAST News

22.10.2025
Updates from EUCAST AFST,
October 2025

21.10.2025

Organization

Public consultations

EUCAST News

Definitions of S, I and R

Clinical breakpoints and dosing

Rapid AST in blood cultures

Expert rules and expected phenotypes

Resistance mechanisms

Guidance documents

SOP

MIC and zone distributions and ECOFFs

AST of bacteria

AST of mycobacteria

AST of fungi

AST of phages

AST of veterinary pathogens

Frequently Asked Questions (FAQ)

Meetings

Rationale documents and publications

Presentations and statistics

launch.

We're getting a fresh new look! Our ne





The European Committee on Antimicrobial
Susceptibility Testing - EUCAST

EUCAST is a standing committee jointly organized by ESCMID, ECDC and European
national bre
Phillips (199
Christian Gi
were Derek
Wootton (20
The EUCAST
Rafael Cant
EUCAST ec
bacteria Gu
EDL for funi
The Genera
China (2025
(next period
EUCAST pr

Κλινικά όρια ευαισθησίας
Δοσολογικά σχήματα
Επιδημιολογικά cut-offs (ECOFs)
Expert rules / Ενδογενής αντοχή
Ανίχνευση μηχανισμών αντοχής
Rational documents (breakpoint background)

Μέθοδος διάχυσης των δίσκων Kirby-Bauer: ερμηνεία

S, I, R



Enterobacteriales^{*} Expert Rules and Expected Phenotypes

Guidance documents

EUCAST Clinical Breakpoint Tables v. 15.0, valid from 2025-01-01 For abbreviations and explanations of breakpoints, see the Notes sheet

Carbapenems ¹	MIC breakpoints (mg/L)			Disk content (μg)	Zone diameter breakpoints (mm)			Notes
	S ≤	R >	ATU		S ≥	R <	ATU	
Doripenem	1	2		10	24	21		1. Some isolates that produce carbapenemase are categorised as susceptible with the current breakpoints and should be reported as tested, i.e. the presence or absence of a carbapenemase does not in itself influence the categorisation of susceptibility. Carbapenemase detection and characterisation are recommended for public health and infection control purposes. For carbapenemase screening, a meropenem screening cut-off of >0.125 mg/L (zone diameter <28 mm) is recommended.
Ertapenem	0.5	0.5		10	23	23		
Imipenem, Enterobacteriales except Morganellaceae	2	4		10	22	19		
Imipenem ² , Morganellaceae	0.001	4		10	50	19		
Imipenem-relebactam, Enterobacteriales except Morganellaceae	2 ³	2 ³		10-25	22	22	20-22	2. The intrinsically low activity of imipenem against <i>Morganella morganii</i> , <i>Proteus</i> spp. and <i>Providencia</i> spp. requires the high exposure of imipenem.
Meropenem (indications other than meningitis)	2	8		10	22	16		3. For susceptibility testing purposes, the concentration of relebactam is fixed at 4 mg/L.
Meropenem (meningitis)	2	2		10	22	22		4. For susceptibility testing purposes, the concentration of vaborbactam is fixed at 8 mg/L.
Meropenem-vaborbactam	8 ⁴	8 ⁴		20-10	20	20	15-19 ^A	A. For isolates in the ATU, if resistant to meropenem report resistant to meropenem-vaborbactam. If not resistant to meropenem, investigate further.

Monobactams	MIC breakpoints (mg/L)			Disk content (μg)	Zone diameter breakpoints (mm)			Notes
	S ≤	R >	ATU		S ≥	R <	ATU	
Aztreonam ¹	1	4		30	26	21		1. The aztreonam breakpoints for <i>Enterobacteriales</i> will detect clinically important resistance mechanisms (including ESBL). Some isolates that produce beta-lactamases are susceptible to aztreonam with these breakpoints and should be reported as tested, i.e. the presence or absence of an ESBL does not in itself influence the categorisation of susceptibility. ESBL detection and characterisation are recommended for public health and infection control purposes.
Aztreonam-avibactam	4 ²	4 ²		30-20	26	26	22-24	

Fluoroquinolones	MIC breakpoints (mg/L)			Disk content (μg)	Zone diameter breakpoints (mm)			Notes
	S ≤	R >	ATU		S ≥	R <	ATU	
Ciprofloxacin, <i>Salmonella</i> spp. ¹	0.06	0.06			Note ^A	Note ^A		1. There is clinical evidence for ciprofloxacin to indicate a poor response in systemic infections caused by <i>Salmonella</i> spp. with any detectable fluoroquinolone resistance mechanisms. The available data relate mainly to <i>Salmonella</i> Typhi but there are also case reports of poor response with other <i>Salmonella</i> species.
Ciprofloxacin (indications other than meningitis)	0.25	0.5	0.5	5	26	22	22-24	
Ciprofloxacin (meningitis) ²	0.125	0.125			Note ^B	Note ^B		
Pefloxacin (screen only)	NA	NA		5	24 ^{A,B,C}	24 ^{A,B,C}		
Delafloxacin, <i>E. coli</i>	0.125	0.125			Note ^D	Note ^D		2/B. In meningitis, where all fluoroquinolone resistance mechanisms must be excluded, either perform an MIC test, or infer susceptibility from the pefloxacin 5 μg screening test. 3. Fluoroquinolone breakpoints are available for other agents. A. Tests with a ciprofloxacin 5 μg disk will not reliably exclude all fluoroquinolone resistance mechanisms in <i>Salmonella</i> spp. Perform an MIC test, or infer susceptibility from the pefloxacin 5 μg screening test. C. The pefloxacin screening test can also be used to detect fluoroquinolone resistance mechanisms in other <i>Enterobacteriales</i> such as <i>E. coli</i> , <i>K. pneumoniae</i> and <i>Shigella</i> spp. D. A disk diffusion test awaits action from the responsible pharmaceutical company.
Levofloxacin	0.5	1		5	23	19		
Moxifloxacin, Enterobacteriales except Morganella morganii, Proteus spp. and Serratia spp. ³	0.25	0.25		5	22	22		
Nalidixic acid (screen only)	NA	NA			NA	NA		
Norfloxacin (uncomplicated UTI only)	0.5	0.5		10	24	24		
Ofloxacin	0.25	0.5		5	24	22		

The European Committee on Antimicrobial Susceptibility Testing. Breakpoint tables for interpretation of MICs and zone diameters. Version 15.0, 2025. <http://www.eucast.org>.

Methodology - EUCAST rapid antimicrobial susceptibility testing (RAST)
directly from positive blood culture bottles.

Version 6.0

March 2025

EUCAST RAST Breakpoint Tables version 8.0 (2025-03-18)

European Committee on Antimicrobial Susceptibility Testing
**Zone diameter breakpoint tables for rapid antimicrobial susceptibility testing (RAST)
directly from blood culture bottles**
Version 8.0, valid from 2025-03-18

Screening for ESBL and carbapenemases in *Escherichia coli*, *Klebsiella pneumoniae* and *Salmonella enterica* for epidemiological purposes as part of the RAST procedure.

EUCAST Guidelines for detection of resistance mechanisms and specific resistance of clinical and/or epidemiological importance using EUCAST rapid antimicrobial susceptibility testing (RAST) directly from positive blood culture bottles.

Version 3.0

March 2025

EUCAST, Rapid AST

- Χρήση με φιάλες BACTEC (BD), BacT/ALERT (bioMérieux) και VersaTREK (Thermo Fisher)
- 0 – 18 h μετά τη θετικοποίηση της φιάλης
- Απευθείας ενοφθαλμισμός (MH, MH-F) με 125±25 μl από τη θετική φιάλη
- Το μικροβιακό είδος **πρέπει να είναι γνωστό** πριν την ερμηνεία των αποτελεσμάτων

Table 1. Incubation conditions for antimicrobial susceptibility test plates.

Organism	Incubation time	Medium	Incubation
<i>Escherichia coli</i> <i>Klebsiella pneumoniae</i> <i>Salmonella enterica</i> <i>Acinetobacter baumannii</i> <i>Staphylococcus aureus</i> <i>Enterococcus faecalis</i> <i>Enterococcus faecium</i>	4, 6 and 8 hours 16-20 hours	MH	35±1°C in air
<i>Pseudomonas aeruginosa</i>	6 and 8 hours 16-20 hours	MH	35±1°C in air

Ταινίες διαβαθμισμένης συγκέντρωσης αντιβιοτικών



1
Prepare inoculum to
0.5 – 1 McFarland



2
Inoculate plate,
allowing drying time of
15 – 20 minutes



3
Position ETEST
strip and incubate
overnight



4
Read and interpret
results

The MIC value is read from the scale in terms of $\mu\text{g/mL}$ at complete inhibition of bacterial growth, where the pointed end of the ellipse intersects the strip.



SET UP TIME:

< 5 Minutes

TIME TO RESULT:


16 – 20 hours for most organisms

Τιμή MIC



Μέθοδος μακροαραιώσεων σε ζωμό

1. Obtain isolated colonies of bacterial strain to test.




2. Combine 4-5 colonies and culture overnight in rich media broth.





Broth dilution method for measuring minimum inhibitory concentration of antibiotics

3. After overnight incubation shown at left, add rich broth with appropriate dilution series of test antibiotic to test tubes. Example concentrations (mg/L) are shown below. Inoculate bacteria to a final density of 5×10^5 cfu/ml.



No bacteria; broth control

4. Plate aliquot of growth control (i.e., no antibiotic added) to verify cfu/ml counts of viable bacteria. Incubate overnight and count colonies.




5. After overnight incubation, check cultures for growth. The MIC is the lowest concentration of antibiotic that prevents visible growth. In this example, the MIC is 64 mg/L.


Correct inoculum
 5×10^5 CFU/ml

Μέθοδος μικροαραίωσης σε ζωμό (BMD)



1. Preparation of test inoculum



2. Preparation of different dilutions of antimicrobial agent



3. Inoculation on 96 well plate




Kept for incubation at 37 °C for 18 hours

4. Results - Determination of Minimum Inhibitory Concentration (MIC)



(Turbidity of the sample is determined)



The lowest concentration of antimicrobial agent that is capable of inhibiting bacterial growth is called MIC.




EUCAST EUROPEAN COMMITTEE
ON ANTIMICROBIAL
SUSCEPTIBILITY TESTING
European Society of Clinical Microbiology and Infectious Diseases



Version 5.0, January 2024

EUCAST reading guide for broth microdilution

- BMD is the reference method for AST of rapidly growing aerobic bacteria, except for mecillinam and fosfomycin, where agar dilution is the reference method.
- EUCAST recommends testing according to the International Standard ISO 20776-1 (with the use of MH-F broth for fastidious organisms).





Μέθοδος αραίωσης σε άγαρ (agar dilution)



mecillinam and fosfomycin:
agar dilution is
the reference method

Fosfomycin Agar Dilution Panel, Liofilchem



Fosfomycin MIC range: 0.25 - 256 µg/mL

A	256	128	64	32
B	16	8	4	2
C	1	0.5	0.25	Control

Growth-control: No antimicrobial agent in the well.

Αυτοματοποιημένα συστήματα




Detection of AMR facilitated by Advanced Expert System (AES).

The instruments read the kinetics of the growth using an optical system with photometer/multichannel fluorimeter readings to assess fluorescence / turbidity /colorimetric signals



- Each of the systems has inherent advantages and limitations
- Results can vary widely by antimicrobial drugs, software versions, and cards
- Some of the systems are not reliable for correct categorization of susceptibility for certain drugs, leading to wrong classifications.
- Software updates and synchronization of breakpoints according to the current standards are mandatory.
- Panels usually contain only several concentrations of each antimicrobial agent, and the resulting MIC is not always given as an exact value.

Advantages and disadvantages of the common methods of antimicrobial susceptibility testing

Method	Advantage	Disadvantage	Comments
Broth dilution 	Well-standardised Harmonised Commercially available tests are easy to perform	Time-consuming Individual mistakes	Quantitative **
Agar Dilution	Well-standardised Suitable for testing a large number of isolates	Time-consuming Limited concentration of antimicrobial agents	Quantitative Possible automation in part
Disk diffusion 	Simple to perform Low cost Simple and fast interpretation The high number of test antibiotics per test High flexibility in antibiotic selection Detection of resistance patterns Mass use and the possibility of automatisisation A number of a different use (AST, identification, screening, etc.) Detection of heteroresistant population or contamination	Time-consuming No MIC value The inability for some antibiotics to be tested	Qualitative *
Gradient test	Convenient and flexible Simple to perform Does not require expertise Detection of resistance patterns	Relatively expensive Relatively long incubation	Quantitative
Automated systems	Simple to perform	Relatively expensive	Semi-quantitative ***

Όχι τιμές MIC

Antibiotic resistance threats

Carbapenem-resistant
Enterobacterales (CRE)

Extended-spectrum
cephalosporin resistance in
Enterobacterales suggestive of
extended-spectrum β -lactamase
(ESBL) production


Multidrug-resistant (MDR)
Pseudomonas aeruginosa

Methicillin-resistant
Staphylococcus aureus
(MRSA)

Carbapenem-resistant
Acinetobacter species

Vancomycin-resistant
Enterococcus (VRE)

Αλγόριθμος για φαινοτυπική ανίχνευση ESBLs σε *Enterobacterales*



¹If ceftazidime has an MIC >8 mg/L,
perform cefepime +/- clavulanic acid
confirmation test

If confirmation with cefepime +/- clavulanic acid is still
indeterminate, genotypic testing is required.


Φαινοτυπική ανίχνευση ESBL

Bionumber: 0405610450106601
Organism Quantity:

Selected Organism: Escherichia coli

Antimicrobial	MIC	Interpretation	Antimicrobial	MIC	Interpretation
ESBL	POS	+	Ertapenem	≤ 0.12	S
Temocillin			Imipenem	≤ 0.25	S
Urine	≤ 4	I	Meropenem		
Other	≤ 4		Meningitis	≤ 0.25	S
Ampicillin	≥ 32	R	Other	≤ 0.25	S
Amoxicillin/Clavulanic Acid	8	S	Amikacin	2	S
Piperacillin/Tazobactam	≤ 4	S	Gentamicin	≤ 1	S
Cefuroxime	≥ 64	R	Tobramycin	≤ 1	S
Cefuroxime Axetil	≥ 64	R	Ciprofloxacin		
Cefoxitin	≤ 4	IE	Meningitis	0.5	
Cefditoren			Other	0.5	I
Cefixime	≥ 4	R	Levofloxacin	1	I
Cefotaxime			Moxifloxacin	0.5	R
Meningitis	≥ 64	R	Minocycline	(-)	(-)
Other	≥ 64	R	Tetracycline	(-)	(-)
Ceftazidime	8	R	Tigecycline	≤ 0.5	S
Ceftriaxone			Fosfomycin		
Meningitis	≥ 64	R	Oral	≤ 16	
Other	≥ 64	R	Other	≤ 16	S
Ceftazidime/Avibactam	≤ 0.12	S	Nitrofurantoin	≤ 16	S
Ceftolozane/Tazobactam	≤ 0.25	S	Chloramphenicol	4	IE
Cefepime	2	I	Colistin	≤ 0.5	S
Aztreonam	16	R	Trimethoprim/ Sulfamethoxazole	≤ 20	S


DDST



CDT



Etest ESBL





≥ 5 mm increase
in inhibition zone

MIC ratio ≥ 8
or if a
phantom zone
or deformed
ellipse is
present

Φαινοτυπική ανίχνευση ESBL


Etest ESBL



EUCAST guidelines for detection of resistance mechanisms and specific resistances of clinical and/or epidemiological importance. Version 2.01. July 2017

Χρωματομετρικές μέθοδοι

β LACTA™ test



directly with isolated colonies or with bacterial pellets from positive blood cultures or urines

Se: 98%
Sp: 100%


Χρωμογόνα υλικά

Liofilchem®
Chromatic ESBL



Escherichia coli DSM 22311


Brilliance™
ESBL Agar



CHROMID® ESBL



CHROMagar™
ESBL



Del Corpo O et al. Clin Microbiol Infect. 2023:S1198-743X(23)00425-1.

CRE-Vitek reports

Organism Quantity:

Selected Organism: Klebsiella pneumoniae


Antimicrobial	MIC	Interpretation	Antimicrobial	MIC	Interpretation
ESBL	NEG	-	Ertapenem	≥ 8	R
Temocillin			Imipenem	≥ 16	R
Urine	≥ 32	R	Meropenem		
Other	≥ 32		Meningitis	≥ 16	R
Ampicillin	≥ 32	R	Other	≥ 16	R
Amoxicillin/Clavulanic Acid	≥ 32	R	Amikacin	16	R
Piperacillin/Tazobactam	≥ 128	R	Gentamicin	≥ 16	R
Cefuroxime	≥ 64	R	Tobramycin	≥ 16	R
Cefuroxime Axetil	≥ 64	R	Ciprofloxacin		
Cefoxitin	≥ 64	IE	Meningitis	≥ 4	
Cefditoren			Other	≥ 4	R
Cefixime	≥ 4	R	Levofloxacin	≥ 8	R
Cefotaxime			Moxifloxacin	≥ 8	R
Meningitis	≥ 64	R	Minocycline	(-)	(-)
Other	≥ 64	R	Tetracycline	(-)	(-)
Ceftazidime	≥ 64	R	Tigecycline	(2)	-
Ceftriaxone			Fosfomycin		
Meningitis	≥ 64	R	Oral	128	R
Other	≥ 64	R	Other	128	R
Ceftazidime/Avibactam	4	(S)	Nitrofurantoin		
Ceftolozane/Tazobactam	≥ 32	R	Chloramphenicol	≥ 64	IE
Cefepime	≥ 32	R	Colistin	≥ 16	(R)
Aztreonam	≥ 64	R	Trimethoprim/ Sulfamethoxazole	≥ 320	R

Organism Quantity:

Selected Organism: Klebsiella pneumoniae

Antimicrobial	MIC	Interpretation	Antimicrobial	MIC	Interpretation
ESBL	NEG	-	Ertapenem	≥ 8	R
Temocillin			Imipenem	≥ 16	R
Urine	≥ 32	R	Meropenem		
Other	≥ 32		Meningitis	≥ 16	R
Ampicillin	≥ 32	R	Other	≥ 16	R
Amoxicillin/Clavulanic Acid	≥ 32	R	Amikacin	4	(S)
Piperacillin/Tazobactam	≥ 128	R	Gentamicin	≥ 16	R
Cefuroxime	≥ 64	R	Tobramycin	8	R
Cefuroxime Axetil	≥ 64	R	Ciprofloxacin		
Cefoxitin	≥ 64	IE	Meningitis	≥ 4	
Cefditoren			Other	≥ 4	R
Cefixime	≥ 4	R	Levofloxacin	≥ 8	R
Cefotaxime			Moxifloxacin	≥ 8	R
Meningitis	≥ 64	R	Minocycline	(-)	(-)
Other	≥ 64	R	Tetracycline	(-)	(-)
Ceftazidime	≥ 64	R	Tigecycline	(1)	-
Ceftriaxone			Fosfomycin		
Meningitis	≥ 64	R	Oral	≤ 16	
Other	≥ 64	R	Other	≤ 16	(S)
Ceftazidime/Avibactam	≥ 16	R	Nitrofurantoin		
Ceftolozane/Tazobactam	≥ 32	R	Chloramphenicol	16	IE
Cefepime	≥ 32	R	Colistin	≤ 0.5	(S)
Aztreonam	≥ 64	R	Trimethoprim/ Sulfamethoxazole	≤ 20	(S)


Ανίχνευση παραγωγής καρβαπενεμασών στα Εντεροβακτηριακά: πληροφορίες από το αντιβιογράμμα



KPC / ESBP +/- [Caz-avi S]


MBL / ESBP + [Caz-avi R]

OXA-48 / ESBP +



MBL

ESBP negative



OXA-48

ESBP negative

Ανίχνευση μηχανισμών αντοχής: CRE

Clinical breakpoints and **screening cut-off values** for carbapenemase-producing Enterobacteriaceae (according to EUCAST methodology)

Carbapenem	MIC (mg/L)		Disk diffusion zone diameter (mm) with 10 µg disks	
	S/I breakpoint	Screening cut-off	S/I breakpoint	Screening cut-off
Meropenem ¹	≤2	>0.125	≥22	<28 ²
Ertapenem ³	≤0.5	>0.125	≥25	<25

¹Best balance of sensitivity and specificity


²Isolates with **25-27 mm only need to be investigated for carbapenemase-production if they are resistant to piperacillin-tazobactam and/or temocillin** (temocillin contributes more to the specificity).

Investigation for carbapenemases is always warranted if zone diameter of meropenem is <25 mm.

³High sensitivity but low specificity. Can be used as an alternative screening agent, but isolates with ESBL and AmpC may be resistant without having carbapenemases.

Phenotypic methods


- combination disk test methods
- colorimetric assays based on hydrolysis of carbapenems




- other methods detecting carbapenem hydrolysis
- lateral flow assays

Ανίχνευση μηχανισμών αντοχής: CRE

Φαινοτυπικός αλγόριθμος για ανίχνευση παραγωγής καρβαπενεμασών στα Εντεροβακτηριακά




Use of the modified cloverleaf (Hodge) test is **not currently recommended** as results are difficult to interpret, the specificity is poor, and in some cases the sensitivity is also suboptimal



¹Combination of several carbapenemases can also contribute to no synergy – e.g. MBL and KPC in combination.

Molecular testing is usually necessary in such cases

²High-level temocillin resistance (>128 mg/L, zone diameter <11 mm) is a phenotypic marker of OXA-48



Journal of Antimicrobial Chemotherapy (2008) **62**, 1257–1260
doi:10.1093/jac/dkn364
Advance Access publication 4 September 2008

JAC

First occurrence of KPC-2-possessing *Klebsiella pneumoniae* in a Greek hospital and recommendation for detection with boronic acid disc tests

Athanassios Tsakris^{1*}, Ioulia Kristo², Aggeliki Poulou³, Fani Markou³,
Alexandros Ikonomidis² and Spyros Pournaras²

JOURNAL OF CLINICAL MICROBIOLOGY, Feb. 2009, p. 362–367
0095-1137/09/\$08.00+0 doi:10.1128/JCM.01922-08
Copyright © 2009, American Society for Microbiology. All Rights Reserved.

Vol. 47, No. 2

Evaluation of Boronic Acid Disk Tests for Differentiating KPC-Possessing *Klebsiella pneumoniae* Isolates in the Clinical Laboratory[▽]

Athanassios Tsakris,^{1*} Ioulia Kristo,² Aggeliki Poulou,³ Katerina Themeli-Digalaki,⁴
Alexandros Ikonomidis,² Dimitra Petropoulou,⁵ Spyros Pournaras,² and Danai Sofianou⁶

JOURNAL OF CLINICAL MICROBIOLOGY, Nov. 2009, p. 3420–3426
0095-1137/09/\$12.00 doi:10.1128/JCM.01314-09
Copyright © 2009, American Society for Microbiology. All Rights Reserved.

Vol. 47, No. 11

Use of Boronic Acid Disk Tests To Detect Extended-Spectrum β -Lactamases in Clinical Isolates of KPC Carbapenemase-Possessing *Enterobacteriaceae*[▽]

Athanassios Tsakris,^{1*} Aggeliki Poulou,² Katerina Themeli-Digalaki,³ Evangelia Voulgari,¹
Theodore Pittaras,¹ Danai Sofianou,⁴ Spyros Pournaras,⁵ and Dimitra Petropoulou⁶

Journal of
Antimicrobial
Chemotherapy

J Antimicrob Chemother 2010; **65**: 1319–1321
doi:10.1093/jac/dkq124 Advance Access publication 15 April 2010

Inhibitor-based methods for the detection of KPC carbapenemase-producing *Enterobacteriaceae* in clinical practice by using boronic acid compounds

Spyros Pournaras¹, Aggeliki Poulou² and Athanassios Tsakris^{3*}

JOURNAL OF CLINICAL MICROBIOLOGY, Aug. 2011, p. 2804–2809
0095-1137/11/\$12.00 doi:10.1128/JCM.00666-11
Copyright © 2011, American Society for Microbiology. All Rights Reserved.

2011;49(8):2804-9.

Vol. 49, No. 8

Comparative Evaluation of Combined-Disk Tests Using Different Boronic Acid Compounds for Detection of *Klebsiella pneumoniae* Carbapenemase-Producing *Enterobacteriaceae* Clinical Isolates[▽]

Athanassios Tsakris,^{1*} Katerina Themeli-Digalaki,² Aggeliki Poulou,^{1,3} Georgia Vrioni,¹
Evangelia Voulgari,¹ Vasiliki Koumaki,¹ Antonella Agodi,⁴
Spyros Pournaras,⁵ and Danai Sofianou⁶



2013;51(9):2986-90.

A Combined Disk Test for Direct Differentiation of Carbapenemase-Producing *Enterobacteriaceae* in Surveillance Rectal Swabs

Spyros Pournaras,^a Olympia Zarkotou,^{a,b} Aggeliki Poulou,^{a,c} Ioulia Kristo,^d Georgia Vrioni,^a Katerina Themeli-Digalaki,^b
Athanassios Tsakris^a

Department of Microbiology, Medical School, University of Athens, Athens, Greece^a; Department of Microbiology, Tzaneio General Hospital, Piraeus, Greece^b; Department of Microbiology, Serres General Hospital, Serres, Greece^c; Department of Microbiology, Faculty of Medicine, University of Thessaly, Larissa, Greece^d



2014;52(5):1483-9.


Modified CLSI Extended-Spectrum β -Lactamase (ESBL) Confirmatory Test for Phenotypic Detection of ESBLs among *Enterobacteriaceae* Producing Various β -Lactamases

Aggeliki Poulou,^{a,b} Evgenia Grivakou,^a Georgia Vrioni,^a Vassiliki Koumaki,^a Theodoros Pittaras,^a Spyros Pournaras,^a
Athanassios Tsakris^a

Department of Microbiology, Medical School, University of Athens, Athens, Greece^a; Department of Microbiology, General Hospital of Serres, Serres, Greece^b

Φαινοτυπική ανίχνευση καρβαπενεμασών στα Εντεροβακτηριακά: χρήση αναστολέων καρβαπενεμασών

CDT



KPC

CDT






MBL



KPC+MBL

Φαινοτυπική ανίχνευση καρβαπενεμασών στα Εντεροβακτηριακά



VIM (-)	KPC (-)	OXA-48 (+)	CTX-M (+)	TEM (+)	NDM (-)	SHV (+)
---------	---------	-------------------	------------------	----------------	---------	----------------

E. coli

Antimicrobial	MIC	Interpretation	Antimicrobial	MIC	Interpretation
ESBL	POS	+	Ertapenem	0.25	S
Temocillin			Imipenem	≤ 0.25	S
Urine	≥ 32	R	Meropenem		
Other	≥ 32		Meningitis	≤ 0.25	S
Ampicillin	≥ 32	R	Other	≤ 0.25	S
Amoxicillin/Clavulanic Acid	≥ 32	R	Amikacin	2	S
Piperacillin/Tazobactam	≥ 128	R	Gentamicin	≤ 1	S
Cefuroxime	≥ 64	R	Tobramycin	≤ 1	S
Cefuroxime Axetil	≥ 64	R	Ciprofloxacin		
Cefoxitin	≤ 4	IE	Meningitis	0.5	
Cefditoren			Other	0.5	I
Cefixime	≥ 4	R	Levofloxacin	1	I
Cefotaxime			Moxifloxacin	0.5	R
Meningitis	≥ 64	R	Minocycline	(-)	(-)
Other	≥ 64	R	Tetracycline	(-)	(-)
Ceftazidime	32	R	Tigecycline	≤ 0.5	S
Ceftriaxone			Fosfomycin		
Meningitis	≥ 64	R	Oral	≤ 16	
Other	≥ 64	R	Other	≤ 16	S
Ceftazidime/Avibactam	≤ 0.12	S	Nitrofurantoin		
Ceftolozane/Tazobactam	8	R	Chloramphenicol	16	IE
Cefepime	≥ 32	R	Colistin	≤ 0.5	S
Aztreonam	≥ 64	R	Trimethoprim/ Sulfamethoxazole	≥ 320	R



(-)= Susceptibility testing not recommended; species is a poor target for therapy IE= Insufficient Evidence that species is good target for therapy; MIC may be reported without interpretation



Φαινοτυπική ανίχνευση παραγωγής καρβαπενεμασών στα Εντεροβακτηριακά


Χρωματομετρικές μέθοδοι (colorimetric tests)

Carba NP Test



- (1) Revealing solution
(internal negative control)
- (2) Revealing solution + imipenem

(b)




Sensitivity: 73–100% for most carbapenemases
Performed poorly for OXA-48 enzyme

Nordmann P, Poirel L, Dortet L. Emerg Infect Dis. 2012;18:1503-7.
Dortet L, Poirel L, Nordmann P. Antimicrob Agents Chemother. 2012;56:6437-40.
Dortet L et al. J Antimicrob Chemother. 2015 ; 70):3014-22.
Kabir MH et al. J Antimicrob Chemother. 2016 ;71:1213-6.
Noël A et al. J Clin Microbiol. 2017;55(2):510-518
Pasteran F et al. J Clin Microbiol. 2015 ;53(6):1996-8.



Blue-Carba Test

indicator bromothymol blue



Screening test: Carbapenemase Inactivation Methods (mCIM/eCIM)

The Carbapenemase Inactivation Methods (mCIM/eCIM) are Currently Recommended by CLSI



Tamma PD, Simner PJ. J Clin Microbiol. 2018;56(11):e01140–18.
<https://asm.org/Articles/2019/May/Opening-the-Black-Box-Phenotypic-Carbapenemase-Det>


van der Zwaluw K et al. PLoS One. 2015;10(3):e0123690
 Yamada K et al. J Microbiol Methods. 2016;128:48-51.
 Tijet N, Patel SN, Melano RG. J Antimicrob Chemother. 2016 ;71(1):274-6.

Vitek reports

Bionumber: 0201010103500210
Organism Quantity:

Selected Organism: Acinetobacter baumannii complex

Antimicrobial	MIC	Interpretation	Antimicrobial	MIC	Interpretation
ESBL			Imipenem	≥ 16	R
Temocillin			Meropenem		
Ampicillin	(-)	(-)	Meningitis	≥ 16	R
Amoxicillin/Clavulanic Acid			Other	≥ 16	R
Piperacillin/Tazobactam	≥ 128	IE	Amikacin	≥ 64	R
Cefuroxime	(-)	(-)	Gentamicin	≥ 16	R
Cefuroxime Axetil	(-)	(-)	Tobramycin	≥ 16	R
Cefoxitin	(-)	(-)	Ciprofloxacin	≥ 4	R
Cefditoren			Levofloxacin	4	R
Cefixime	(-)	(-)	Moxifloxacin	(-)	(-)
Cefotaxime	(-)	(-)	Minocycline	2	IE
Ceftazidime	(-)	(-)	Tetracycline	(-)	(-)
Ceftriaxone	(-)	(-)	Tigecycline	4	IE
Ceftazidime/Avibactam			Fosfomycin	(-)	(-)
Ceftolozane/Tazobactam			Nitrofurantoin	(-)	(-)
Cefepime	(-)	(-)	Chloramphenicol	(-)	(-)
Aztreonam	(-)	(-)	Colistin	≤ 0.5	S
Ertapenem			Trimethoprim/ Sulfamethoxazole	≥ 320	R



Organism Quantity:

Selected Organism: Pseudomonas aeruginosa

Comments:	Colistin S results / Gram negative bacilli: perform an alternative method of testing prior to reporting results.



Identification Information			
Organism Origin	Technologist		
Selected Organism	Pseudomonas aeruginosa		
	Entered: Sep 15, 2023 13:42 EEST		By: LabTech
Analysis Messages:			

Susceptibility Information	Card:	AST-N318	Lot Number:	7582293203	Expires:	Mar 7, 2024 12:00 EET
	Status:	Final	Analysis Time:	9.27 hours	Completed:	Sep 15, 2023 22:42 EEST

Antimicrobial	MIC	Interpretation	Antimicrobial	MIC	Interpretation
+Ampicillin		(-)	Meropenem		
Ampicillin/Sulbactam			Meningitis	≥ 16	R
Piperacillin/Tazobactam	≥ 128	R	Other	≥ 16	R
Cefoxitin			Amikacin	≥ 64	R
+Cefotaxime		(-)	Gentamicin	≥ 16	IE
Ceftazidime	≥ 64	R	Ciprofloxacin	≥ 4	R
Ceftriaxone	(-)	(-)	Levofloxacin	≥ 8	R
+Ceftolozane/Tazobactam			Tigecycline	(-)	(-)
Cefepime	≥ 64	R	Fosfomycin	(-)	(-)
Aztreonam	≥ 64	R	Colistin	2	S
+Ertapenem		(-)	Trimethoprim/ Sulfamethoxazole	(-)	(-)
Imipenem	≥ 16	R			

P. aeruginosa και αντοχή στις καρβαπενέμες

Multidrug-resistant (MDR)
Pseudomonas aeruginosa



Τάξη	Ένζυμα
A	KPC GES
D	OXA (OXA-40, -50 -198)
B	VIM IMP NDM SPM GIM, AIM, FIM

P. aeruginosa και φαινοτυπική ανίχνευση MBL

Δεν υπάρχουν προτυποποιημένες δοκιμασίες



Οι δοκιμασίες στηρίζονται στη
χρήση χηλικών παραγόντων
EDTA, MPA, DPA

Η απόδοση εξαρτάται από

- μεθοδολογία
(DDST, CDT, Etest MBL)
- β-λακταμικό υπόστρωμα
- είδος του αναστολέα
- ποσότητα του αναστολέα
- χαρακτηριστικά των
στελεχών

**All phenotypic combined disk tests lacked either sensitivity
or specificity for the detection of MBL in *P. aeruginosa***

CDT



EDTA 292 µg (10 µl EDTA 0.1 M)

Δd > 4 mm


EDTA 750 µg, Δd ≥ 6 mm

100% SN και SP

DPA 835 µg, Δd ≥ 5 mm

SN=93,2% SP= 97,1%

DDST




EDTA 292 µg

5 µl MPA 1.4 Mm

DPA 250 µg


Etest MBL, bioMerieux



IP/ IPI ≥ 8

παρουσία ζώνης-φάντασμα

παραμόρφωση των ελλείψεων IP ή IPI




For *P. aeruginosa* the MBL Etest® as well as disk-based assays have been used for several decades, but are hampered **by poor specificity.**

EUCAST guidelines for detection of resistance mechanisms and specific resistances of clinical and/or epidemiological importance

Review

The Revolution of Lateral Flow Assay in the Field of AMR Detection


Hervé Boutal ¹, Christian Moguet ¹, Lilas Pommès ¹, Stéphanie Simon ¹, Thierry Naas ^{2,3,4} and Hervé Volland ^{1,*}



LFIA performance relies mostly on antibody affinity and specificity



an ideal PoC test




PBC

Se: 98%


Sp: 100%

Del Corpo O et al. Clin Microbiol Infect. 2023 Sep 16:S1198-743X(23)00425-1.

Εντεροβακτηριακά με μειωμένη ευαισθησία στις καρβαπενέμες



Lateral flow assay formats








Targets	Technology	Product	Other identifications
Beta Lactamases	Electrochemistry	BL-RED 25	
OXA-48-like	Lateral flow	OXA-48 K-SeT	
	Lateral flow	O.K.N.V.I. RESIST-5	KPC, NDM, VIM, IMP
	Lateral flow	RESIST-3 O.O.K. K-SeT	OXA-163-like, KPC
	Lateral flow	RESIST-3 O.K.N. K-SeT	KPC, NDM
	Reagents	RESIST-BC	Carbapenemases: KPC, NDM, VIM, IMP, OXA-40/58
	Reagents	ReSCape	Carbapenemases: KPC, NDM, VIM, IMP, OXA-40/58
KPC	Lateral flow	O.K.N.V.I. RESIST-5	OXA-48-like, NDM, VIM, IMP
	Lateral flow	RESIST-3 O.O.K. K-SeT	OXA-48-like, OXA-163-like
	Lateral flow	RESIST-3 O.K.N. K-SeT	OXA-48-like, NDM
	Reagents	RESIST-BC	Carbapenemases: OXA-48-like, NDM, VIM, IMP, OXA-40/58
	Reagents	ReSCape	Carbapenemases: OXA-48-like, NDM, VIM, IMP, OXA-40/58
	Lateral flow	RESIST ACINETO	OXA-23-like, OXA-40/58
OXA-163-like	Lateral flow	RESIST-3 O.O.K. K-SeT	OXA-48-like, KPC
NDM	Lateral flow	O.K.N.V.I. RESIST-5	OXA-48-like, KPC, VIM, IMP
	Lateral flow	RESIST-3 O.K.N. K-SeT	OXA-48-like, KPC
	Reagents	RESIST-BC	Carbapenemases: OXA-48-like, KPC, VIM, IMP, OXA-40/58
	Reagents	ReSCape	Carbapenemases: OXA-48-like, KPC, VIM, IMP, OXA-40/58
	Lateral flow	RESIST ACINETO	OXA-23-like, OXA-40/58
	Lateral flow	OXA-23 K-SeT	OXA-40/58, NDM
VIM	Lateral flow	O.K.N.V.I. RESIST-5	OXA-48-like, KPC, NDM, IMP
	Reagents	RESIST-BC	Carbapenemases: OXA-48-like, KPC, NDM, IMP, OXA-40/58
	Reagents	ReSCape	Carbapenemases: OXA-48-like, KPC, NDM, IMP, OXA-40/58
IMP	Lateral flow	O.K.N.V.I. RESIST-5	OXA-48-like, KPC, NDM, VIM
	Lateral flow	IMP K-SeT	
	Reagents	RESIST-BC	Carbapenemases: OXA-48-like, KPC, NDM, VIM, OXA-40/58
	Reagents	ReSCape	Carbapenemases: OXA-48-like, KPC, NDM, VIM, OXA-40/58
Carbapenemases	Reagents	RESIST-BC	Carbapenemases: OXA-48-like, KPC, NDM, VIM, IMP, OXA-40/58
	Reagents	ReSCape	Carbapenemases: OXA-48-like, KPC, NDM, VIM, IMP, OXA-40/58
OXA-40/58	Reagents	RESIST-BC	Carbapenemases: OXA-48-like, KPC, NDM, VIM, IMP
	Reagents	ReSCape	Carbapenemases: OXA-48-like, KPC, NDM, VIM, IMP
	Lateral flow	RESIST ACINETO	OXA-23-like, NDM
Group 1 (CTX-M-15,...)	Lateral flow	RESIST CTX-M	Group 9 (CTX-M-14,...)
Group 9 (CTX-M-14,...)	Lateral flow	RESIST CTX-M	Group 1 (CTX-M-15,...)

ACCESSORIES TO ACCELERATE RESIST RANGE RESULTS

PRODUCT NAME	DESCRIPTION
RESIST-BC	Reagents kit for blood culture 1 RBCL solution, 1 MS buffer, 1 Washing buffer
ReSCape	Culture reagent kit on rectal swab 20 CPROBE MEDIUM tubes, 2 Selective Mix vials, 2 Water vials



Se: 98 -100%
Sp: 100%



- Wareham DW, et al. J Clin Microbiol. 2016;54:471-3.
 Pasteran F, et al. J Clin Microbiol. 2016;54:2832-6.
 Glupczynski Y et al. J Antimicrob Chemother. 2017;72(7):1955-1960.
 Greissl C et al. Eur J Clin Microbiol Infect Dis. 2019;38(2):331-335.
 Song W et al. Ann Lab Med. 2020;40(3):259-263.
 Sadek M et al. Diagn Microbiol Infect Dis. 2022;104(4):115761.
 Bouvier M et al. Diagn Microbiol Infect Dis. 2023;107(3):116043.


Ανίχνευση καρβαπενεμασών σε στελέχη *Acinetobacter*

N = 174 well-characterized carbapenemase and non-carbapenemase producing isolates belonging to various *Acinetobacter* species

149 *A. baumannii*, 11 *A. pittii*, four *A. radioresistens*, three *A. ursingii*, two *A. calcoaceticus*, two *A. nosocomialis*, one *A. bereziniae*, one *A. lwoffii* and one *A. junii*.



Negative test
A. baumannii



OXA-40/58 +
A. baumannii



OXA-23 +
A. baumannii



NDM +
A. baumannii



OXA-40/58 +
OXA-23 +
A. baumannii



OXA-40/58 +
NDM +
A. baumannii




OXA-40/58 +
A. pittii




OXA-23 +
A. radioresistens

Overall sensitivity and specificity: 100% and 96%
PPA and NPA 91% and 100%.

After exclusion of *A. pittii* and *A. radioresistens*
results, sensitivity, specificity, PPA and NPA were
100%.



IDENTIFICATION PROCESS FROM DIRECT BLOOD CULTURE: 20–40 min protocol



Bodendoerfer E et al. J Antimicrob Chemother. 2019;74(6):1749-1751.

NG•TEST®/CARBA-5



Boutal H et al. J Antimicrob Chemother. 2018;73:909-15.
Han R et al. Front Microbiol. 2021;11:609856.
Volland H, et al. Antimicrob Agents Chemother. 2019;64(1):e01940-19.

NG•TEST®/CTX-M Multi


Rapid detection of Extended Spectrum Beta-Lactamase (ESBL)



Bernabeu S et al. Diagnostics (Basel). 2020;10(10):764.
Bianco G, et al. J Hosp Infect. 2020;105(2):341-343.
Boattini M et al. Eur J Clin Microbiol Infect Dis. 2021;40(7):1495-1501.

NG•TEST®/MCR-1

Rapid detection of Colistin resistance



Volland H et al. J Clin Microbiol. 2019;57:e01454-18.
Fenwick AJ et al. J Clin Microbiol. 2020;58(4):e01823-19.

LFIA και ανίχνευση φορέων απευθείας από το δείγμα

J Antimicrob Chemother. 2019 Feb 1;74(2):357-359. doi: 10.1093/jac/dky429.

Lateral flow immunochromatographic assay for rapid screening of faecal carriage of carbapenemase-producing Enterobacteriaceae

Charlotte Fauconnier¹, Magali Dodemont², Angélique Depouhon¹, Ahalieyah Anantharajah¹, Alexia Verroken¹, Hector Rodriguez-Villalobos¹

Residual rectal swabs ($n=149$), routinely screened for CPE through selective culture and confirmed by PCR

Protocol: 2.5h incubation of the swab in an enrichment medium containing meropenem followed by **O.K.N. K-SeT** testing after centrifugation.

Strains: 82 OXA-48, 22 NDM, 5 KPC, 1 IMP, 7 VIM, 3 OXA-23

Detected: 76 OXA-48, 19 NDM, 5 KPC

Sensitivity: 96%

Specificity: 100%

LoD: 10^4 - 10^5 cfu/ml

Microorganisms. 2021 Apr 27;9(5):942. doi: 10.3390/microorganisms9050942.

Detection of KPC, NDM and VIM-Producing Organisms Directly from Rectal Swabs by a Multiplex Lateral Flow Immunoassay

Alexandra Vasilakopoulou¹, Polyxeni Karakosta¹, Sophia Vourli¹, Eleni Kalogeropoulou¹, Spyros Pournaras¹

Fecal swabs ($n=30$)

- Culture
- PCR
- **NG-Test® CARBA-5**

Sensitivity:

80% (direct testing)
88% (1h incubation)
92% (2h incubation)
100% (3h incubation)

Specificity: 100%

Each sample was tested by CARBSA-5 five times:


- Direct testing
- Incubation in TSB for 1h, 2h, 3h, 4h

TSB enriched with 0.25mg/l meropenem

Concordant results from LFIA and PCR: 11 KPC, 8 VIM, 6 NDM. [In five cases, co-carriage of KPC and VIM]

Χρωμογόνα υλικά

- Purpose: more rapid detection and identification of MDRO
- The target organisms are characterized by specific enzyme systems that metabolize the substrates to release the chromogen.
- The chromogen can then be visually detected by direct observation of a distinct colour change in the medium (coloured colonies).
- Compared with the use of conventional culture media, the use of chromogenic agar often reduces the costs and labor time.
- Their primary use is for screening of patients colonized with various pathogens.
- The sensitivity and specificity depend on the manufacturer and the type of microorganism detected
- Additional confirmation of the resistant bacteria is sometimes needed.



Advantage

- Mass use and the possibility of automatisation
- Simple to perform
- Simple and fast interpretation

Disadvantage

- Not completely susceptible and specific
- Time-consuming
- Limited spectra or single antibiotic
- Relatively expensive
- Screening only or required confirmatory identification
- No MIC value

Comments

Qualitative with no interpretation criteria (S, I, R)

Χρωμογόνα εκλεκτικά υλικά για CRE


CHROMagar KPC, Colorex KPC (CHROMagar)



Klebsiella, Enterobacter, Citrobacter CR
→ **Metallic blue**
Pseudomonas CR → Cream, translucent
E.coli CR → **dark pink to reddish**

- Μειονέκτημα: έλλειψη ευαισθησίας
- Δυσκολία ανάπτυξης στελεχών με $MIC \leq 4 \mu g/ml$
- Προβληματική η ανίχνευση στελεχών MBL και OXA-48

ChromID Carba (bioMérieux)



- SN: 92,4%, SP: 96, 9%
- Ικανότητα ανίχνευσης στελεχών με χαμηλές MIC¹
- Υψηλότερη ικανότητα ανίχνευσης για KPC-στελέχη vs Brilliance CRE /Colorex KPC (σε υψηλό και σε χαμηλό inoculum)²
- Υψηλότερη ικανότητα ανίχνευσης NDM-στελεχών³


1. Vroni G, et al. J Clin Microbiol. 2012; 50(6):1841-6.
2. Wilkinson KM, et al. J Clin Microbiol. 2012; 50(9):3102-4.
3. Perry JD, et al. J Antimicrob Chemother. 2011; 66(10):2288-94.

ChromID OXA-48 (bioMérieux)



Incubation : 18 hours
K. pneumoniae OXA-48


CHROMID® CARBA SMART Agar




Χρωμογόνα εκλεκτικά υλικά για CRE

Brilliance™ CRE (Oxoid)


- Χαμηλότερη SN έναντι του ID Carba για την ανίχνευση NDM-στελεχών
- Υπεροχή έναντι του CHROMagar KPC (43%)
- Παρόμοια SP



Brilliance™ ESBL Agar/
Brilliance™ CRE Agar



Liofilchem® Chromatic CRE



Liofilchem® Chromatic OXA-48



OXA-48 positive *Enterobacter cloacae*


Έλεγχος φορείας για CRE: καλλιεργητικές μέθοδοι

Home-made τρυβλία MacConkey agar:

- Με ενσωματωμένη καρβαπενέμη
 - ικανοποιητική ευαισθησία
- Με δίσκους καρβαπενεμών:
 - SN/SP εξαρτώνται από το όριο
- Με δίσκους + αναστολείς:
 - **ανίχνευση + διαφοροποίηση**

Χρωμογόνα υλικά

- **chromID CARBA / chromID OXA-48, bioMérieux**
 - Υψηλή ευαισθησία και ειδικότητα
- **Brilliance CRE, Oxoid**
 - Ικανοποιητική ευαισθησία και ειδικότητα
- **CHROMagar KPC/ Colorex KPC, CHROMagar**
 - Ανιχνεύει στελέχη με υψηλού επιπέδου αντοχή





Pournaras S, et al. J Clin Microbiol. 2013; 51(9):2986-90.


Η επίδοση της κάθε μεθόδου εξαρτάται:

- επίπεδο αντοχής στις καρβαπενέμες
- υδρόλυση άλλων β-λακταμών (φάσμα)
- παρουσία άλλων μηχανισμών αντοχής
- γεωγραφική περιοχή, βακτηριακό είδος, δείγμα, χρόνο επώασης, μέθοδο αναφοράς, ορισμό αληθώς θετικών, ...

Σύγκριση μελετών-μεθόδων δύσκολη




mec gene complex *ccr* gene complex



Φαινοτυπικές μέθοδοι

Disk diffusion



CEFOXITIN: a very **sensitive and specific marker** of *mecA/mecC*-mediated methicillin resistance including in heterogeneous expressing strains and **is the agent of choice**.

Broth microdilution




Organism Quantity:		Selected Organism: <u>Staphylococcus aureus</u>			
Comments:	S.aureus - Mupirocin: ECOFF				
Identification Information					
Organism Origin		Technologist			
Selected Organism		Staphylococcus aureus			
		Entered: Oct 7, 2023 17:56 EEST		By: LabTech	
Analysis Messages:					
Susceptibility Information	Card:	AST-P659	Lot Number:	8092252103	Expires: Jan 26, 2024 12:00 EET
	Status:	Final	Analysis Time:	8.48 hours	Completed: Oct 8, 2023 02:22 EEST
Antimicrobial	MIC	Interpretation	Antimicrobial	MIC	Interpretation
Cefoxitin Screen	POS	+	Linezolid	2	S
Benzylpenicillin	>= 0.5	R	Daptomycin	0.25	S
Oxacillin	>= 4	R	Teicoplanin	<= 0.5	S
Ceftaroline			Vancomycin	<= 0.5	S
Pneumonia	0.25	S	Tetracycline	<= 1	S
Other	0.25	S	Tigecycline	<= 0.12	S
Gentamicin	<= 0.5	S	Fusidic Acid	8	R
Levofloxacin	>= 8	R	Mupirocin	<= 1	S
Inducible Clindamycin Resistance	NEG	-	Rifampicin	<= 0.03	S
Erythromycin	0.5	S	Trimethoprim/ Sulfamethoxazole	<= 10	S
Clindamycin	0.25	S			

Detection of PBP2a with latex agglutination kits




*PBP2c is not detected by the majority of commercial assays

Chromogenic media



Χρωμογόνα εκλεκτικά υλικά για MRSA

CHROMID® MRSA
SMART, bioMerieux



ChromID *S. aureus*/
ChromID MRSA, bioMerieux



Brilliance™ MRSA 2
Agar, Oxoid



Brilliance™ MRSA 2 Agar/
Brilliance™ Staph 24 Agar



MRSA Select,
Bio-Rad



ChromAgar MRSA



Variable	CHROMagar [43]	BBL-CHROMagar [43]	CHROMagar MRSA [35]
Sensitivity (%)	82–93	83–94	96–100
Specificity (%)	97–99	98–99	95–97
Turnaround time (h)	24–48	24–48	24–48
Storage for incubation	2–8 °C; dark	2–8 °C; dark	2–8 °C; dark

Variable	MRSASelect [35]	Brilliance MRSA agar [43]	ChromID [43]
Sensitivity (%)	81–93	90–96	83–94
Specificity (%)	92–97	69–87	90–96
Turnaround time (h)	18–24	18–24	
Storage for incubation	Room temperature	2–10 °C; dark	Room temperature

Harbarth S et al. Int J Antimicrob Agents. 2011;37(2):110-7.

Nonhoff C et al. Eur J Clin Microbiol Infect Dis. 2009 Apr;28(4):363-9.

Malhotra-Kumar S et al. J Clin Microbiol. 2010 Apr;48(4):1040-6.

Φαινοτυπικές μέθοδοι για ανίχνευση του MRSA

Phenotypic Method Used	Principle of the Method	On Culture/on Clinical Sample	TAT ¹	Brief Advantages/Disadvantages +/-	No. and Type of <i>S. aureus</i> Analyzed ²	Major Diagnostic Performance ³	Reference
Chromogenic media	Selective media, colorimetric colony detection	Nasal swabs, positive blood culture	24 h	+ Easy to perform and interpret, low-cost – No MSSA identification	72 MRSA 28 MSSA	Sensitivity 100% Specificity 78.6%	[28]
					19 MRSA 29 MSSA	Sensitivity 78.9% Specificity 41.3%	[29]
Rapid BMD	Colorimetric assay, MTT dye as oxidation-reactions indicator	Culture	7 h	+ Easy to perform, rapid, low-cost – Larger validations needed	21 MRSA 19 MSSA	Sensitivity 100% Specificity 94.7%	[30]
Disk diffusion on Mueller-Hinton rapid agar MHR-SIR	Disk diffusion with shorter incubation	Positive blood culture	6–8 h	+ Easy to perform, rapid, low-cost – Larger validations needed	23 <i>S. aureus</i> tested for several antibiotics	CA 97.8% mEs 1.9% VME 0.3%	[31]
					29 MRSA 50 MSSA	TAT 17h shorter, better therapy de-escalation	[32]
Rapid AST	Disk diffusion with shorter incubation	Positive blood culture	4–6 h	+ Easy to perform, rapid, low-cost – Low growth, diameters margins less demarcated	24 MRSA 313 MSSA	Sensitivity 100% Specificity 100%	[33]
					14 MRSA 183 MSSA	Sensitivity 100% Specificity 99.4% ATU 0.5%	[34]
					9 MRSA 212 MSSA	Sensitivity 100% Specificity 99.1% ATU 0.4% ME 0.5%	[35]

Φαινοτυπικές μέθοδοι για ανίχνευση του MRSA

Phenotypic Method Used	Principle of the Method	On Culture/on Clinical Sample	TAT ¹	Brief Advantages/Disadvantages +/–	No. and Type of <i>S. aureus</i> Analyzed ²	Major Diagnostic Performance ³	Reference
VITEK 2	Turbidimetric monitoring of bacterial growth	Culture	4–11 h	+ Automated, rapid, standardized interpretation – Fixed range of testable antibiotics, less accurate for linezolid/vancomycin	72 MRSA 28 MSSA	Sensitivity 97.2% Specificity 100%	[28]
					27 MRSA resistant to linezolid	MEs 85.2%	[36]
					28 MRSA with a vancomycin MIC ≥ 2	MEs 46.4%	[37]
BD Phoenix	Turbidimetric and colorimetric growth detection	Culture	6–16 h	+ Automated, rapid, standardized interpretation – Less accurate for SXT	642 <i>S. aureus</i> analyzed for SXT resistance	CA 91.9% MEs 7.9% mEs 0.3%	[38]
Copan WASP Colibri	Automated disk diffusion	Culture	16 h	+ Automated, flexible, detect heteroresistance and mixed cultures – Should be compared with standard disk diffusion	107 <i>S. aureus</i> tested for several antibiotics	CA 99.9%	[39]
Accelerate PhenoTest™ BC	Single-cell analysis, fluorescent in situ hybridization	Positive blood culture	7 h	+ Automated, monomicrobial call, rapid, identification and AST in the same platform	98 MRSA 86 MSSA	Sensitivity 100% Specificity 98.8% CA 99.5%	[17]
					22 MRSA 2 MSSA	Sensitivity 100% Specificity 100%	[40]
PBP2a latex agglutination assay	Particles with monoclonal antibodies, agglutination reaction	Culture	5 m	+ Rapid, easy to perform and interpret – no PBP2c detection	95 MRSA 10 MSSA	Sensitivity 98.95% Specificity 77.8%	[41]

Sanchini A. Diagnostics (Basel). 2022;12(1):208.

Μοριακές μέθοδοι για ανίχνευση του MRSA απευθείας από το δείγμα: παραδείγματα

Molecular Method Used	Principle of the Method	On Culture/ on Clinical Sample	TAT ¹
Xpert [®] SA Nasal Complete (Cepheid)	Real-time PCR for <i>mecA/C</i> , <i>spa</i> and <i>SCCmec-orfX</i>	Nasal samples	1 h
Xpert [®] MRSA/SA BC Assay (Cepheid)	Real-time PCR for <i>mecA/C</i> , <i>spa</i> and <i>SCCmec-orfX</i>	Positive blood cultures	1.7 h
Xpert [®] MRSA/SA SSTI (Cepheid)	Real-time PCR for <i>mecA/C</i> , <i>spa</i> and <i>SCCmec-orfX</i>	BAL samples	68 m



Molecular Method Used	Principle of the Method	On Culture/ on Clinical Sample	TAT ¹	Brief Advantages/Disadvantages +/-
Hologic Panther Fusion [®] MRSA	PCR and Invader chemistries for <i>mecA/C</i> , <i>gap</i> and <i>SCCmec-orfX</i>	Nasal samples	<3 h	+ Can analyze 350 samples in 8 h – Need comparison with a similar method
MRSA/SA ELITE MGB assay (ELITech-Group)	Real time PCR for <i>mecA/C</i> and a <i>S. aureus</i> specific sequence	Sputum, tracheal aspirate, BAL	<3 h	+ Accurate – Do not target <i>SCCmec-orfX</i> junction
Eazyplex [®] MRSA	LAMP targeting <i>S. aureus</i> , <i>S. epidermidis</i> , <i>mecA/C</i>	Positive blood cultures	1 h	+ Portable; faster than similar methods – Need optimization for CONS





Sanchini A. Diagnostics (Basel). 2022;12(1):208.

VRE, μηχανισμός αντοχής: τροποποίηση στόχου δράσης

Σύνθεση πεπτιδογλυκάνης που φέρει
πλευρικές πενταπεπτιδικές
αλυσίδες που καταλήγουν σε

D-ala-D-lactate
(VanA, VanB, VanD)

D-ala-D-serine
(VanC, VanE,
VanG, VanL, VanN)




Glycopeptide resistance due to *van*-type gene clusters

Resistance level	Acquired							Intrinsic
	High	Variable	Moderate	Low			Low	
Type	VanA	VanB	VanD	VanE	VanG	VanL	VanN	VanC
MIC in mg/L:								
Vancomycin	≥16	≥4	≥64	6–32	12–16	8	8	2–32
Teicoplanin	>8	0,5–1	4–64	0,5	0,5			0,5–1
Expression	Inducible	Inducible	Constitutive/ Inducible (<i>vanD2</i>)	Inducible/ (Constitutive)	Inducible			Constitutive/ Inducible
<i>van</i> ligase gene	<i>vanA</i>	<i>vanB1-B3</i>	<i>vanD1-5</i>	<i>vanE</i>	<i>vanG1-2</i>	<i>vanL</i>	<i>vanN</i>	<i>vanC1-C3</i>
Modified target	D-alanine-D-lactate	D-alanine-D-lactate	D-alanine-D-lactate	D-alanine-D-serine	D-alanine-D-serine	D-alanine-D-serine	D-alanine-D-serine?	D-alanine-D-serine
Conjugative transfer	Yes	Yes	No	No	Yes	No	Yes	No
Location	Plasmid/chromosome on transposon(s)	Plasmid/chromosome ± transposon/ICE ^a	Chromosome	Chromosome	Chromosome on possible ICE	Chromosome?	Plasmid	Chromosome
Distribution	<i>Enterococcus faecium</i> <i>Enterococcus faecalis</i> <i>Enterococcus avium</i> <i>Enterococcus casseliflavus</i> <i>Enterococcus durans</i> <i>Enterococcus gallinarum</i> <i>Enterococcus hirae</i> <i>Enterococcus mundtii</i> <i>Enterococcus raffinosus</i> <i>Staphylococcus aureus</i> <i>Bacillus circulans</i> <i>Oerskovia turbata</i> <i>Arcanobacterium haemolyticum</i> <i>Paenibacillus</i> <i>Rhodococcus</i>	<i>Enterococcus faecium</i> <i>Enterococcus faecalis</i> <i>Enterococcus casseliflavus</i> <i>Enterococcus durans</i> <i>Enterococcus gallinarum</i> <i>Enterococcus hirae</i> <i>Staphylococcus epidermidis</i> <i>Streptococcus</i> <i>Clostridium</i> <i>Ruminococcus</i> <i>Eggerthella</i>	<i>Enterococcus faecium</i> <i>Enterococcus faecalis</i> <i>Enterococcus avium</i> <i>Enterococcus gallinarum</i> <i>Enterococcus raffinosus</i> Non-enterococcal faecal flora	<i>Enterococcus faecalis</i>	<i>Enterococcus faecalis</i> Non-enterococcal faecal flora	<i>Enterococcus faecalis</i>	<i>Enterococcus faecium</i>	<i>Enterococcus gallinarum</i> – <i>vanC1</i> <i>Enterococcus casseliflavus</i> – <i>vanC2/3</i>

Hegstad K et al. Clin Microbiol Infect. 2010;16(6):541-54.

Ανίχνευση VRE: παραδείγματα

Μέθοδος διάχυσης των δίσκων




Isolates must not be reported susceptible before 24 h incubation

Ταινίες διαβαθμισμένης συγκέντρωσης αντιβιοτικών



sharp zone edges

Confirmatory testing with PCR or report **resistant** even if the zone diameter is ≥ 12 mm



Ανίχνευση VRE: παραδείγματα

Αυτοματοποιημένα συστήματα

Bionumber: 532047265773771
Organism Quantity:

Selected Organism: Enterococcus faecium

Antimicrobial	MIC	Interpretation	Antimicrobial	MIC	Interpretation
Cefoxitin Screen			Moxifloxacin	(-)	(-)
Benzylpenicillin	(-)	(-)	Inducible Clindamycin Resistance		
Ampicillin	≥ 32	R	Erythromycin	(-)	(-)
Oxacillin			Clindamycin	(-)	(-)
Imipenem	≥ 16	R	Quinupristin/Dalfopristin	4	R
Gentamicin High Level (synergy)	SYN-S	S	Linezolid	1	S
Streptomycin High Level (synergy)	SYN-R	R	Daptomycin		
Gentamicin			Teicoplanin	≥ 32	R
Ciprofloxacin			Vancomycin	≥ 32	R
Urine	≥ 8	R	Tigecycline	≤ 0.12	S
Other	≥ 8		Trimethoprim/Sulfamethoxazole		

Broth microdilution



Μοριακές μέθοδοι απευθείας από το δείγμα

VanA or VanB resistance in 20 min, from rectal swabs, blood culture or as culture confirmation



VRE Cepheid Xpert® vanA/vanB



VanA SN= 87.1%, SP=99.7%, PPV=98.0%, NPV=97.7%

VanB SN=77.6%



The saving of time to results corresponded to 141 saved isolation days and 292 saved transmission risk days.

Holzkecht BJ et al. New Microbes New Infect. 2017;16:54-59.

Development and validation of a lateral flow immunoassay for rapid detection of VanA-producing enterococci



Saoussen Oueslati^{1,2,†}, Hervé Volland^{3,†}, Vincent Cattoir⁴, Sandrine Bernabeu^{1,2}, Delphine Girlich¹, Duncan Dulac³, Marc Plaisance³, Maxime Laroche⁵, Laurent Dortet^{1,2,6}, Stéphanie Simon³ and Thierry Naas^{1,2,6*}



Article

Rapid Detection of VanA/B-Producing Vancomycin-Resistant Enterococci Using Lateral Flow Immunoassay

Saoussen Oueslati^{1,2,†}, Camille Gonzalez^{1,2,†}, Hervé Volland^{3,*,†}, Vincent Cattoir⁴, Sandrine Bernabeu^{1,2}, Delphine Girlich^{1,2}, Duncan Dulac³, Marc Plaisance³, Laure Boutigny⁵, Laurent Dortet^{1,2}, Stéphanie Simon³ and Thierry Naas^{1,*,†}




They may also be used from positive blood cultures, given a 3.5 h sub-culturing step in the presence of a 5-µg vancomycin disk

Χρωμογόνα εκλεκτικά υλικά για VRE


Brilliance VRE®, Oxoid



Chromatic VRE,
Liofilchem




ChromAgar® VRE



chromID® VRE,
bioMerieux



Bile-esculin-azide agar +
6 mg/L vancomycin



Παράγοντες που επηρεάζουν την ευαισθησία

- χαμηλή πυκνότητα VRE στα κόπρανα
- όταν η ανάγνωση γίνεται στις 24 h (vs 48 h)
- χαμηλές τιμές MIC στη vancomycin (**VanB φαινότυπος**)
 - **SN 25-75% (MIC <16 mg/l) vs 100% for MIC >32 mg/l.**

TABLE 2 Comparison 396 specimens (combined data)

Medium	Performance (% [95% CI]) ^b			
	Sensitivity	Specificity	PPV	NPV
BEAV broth	98.0 (94.7–98.0)	100 (98.9–100)	100 (96.7–100)	99.3 (98.3–99.3)
Enterococcosel (BEAV)	84.8 (80.9–84.8)	100 (98.7–100)	100 (95.3–100)	95.2 (93.9–95.2)
InTray Colorex VRE	91.9 (86.9–94.8)	98.3 (96.6–99.3)	94.8 (89.6–97.9)	97.3 (95.7–98.3)
chromID	94.9 (91.0–95.9)	99.7 (98.3–100)	98.9 (94.8–99.9)	98.3 (97.0–98.7)
VRESelect	91.9 (87.8–92.9)	99.7 (98.3–100)	98.9 (94.4–99.9)	97.4 (96.0–97.7)
HardyCHROM VRE	89.9 (85.6–90.9)	99.7 (98.2–100)	98.9 (94.2–99.9)	96.7 (95.4–97.0)
Spectra VRE	93.9 (89.9–94.9)	99.7 (98.3–100)	98.9 (94.7–99.9)	98.0 (96.7–98.3)

Suwanarat N et al. J Clin Microbiol. 2014;52(11):4039-42.

Wijesuriya TM et al. J Clin Microbiol. 2014 Aug;52(8):2829-33.


Τα χρωμογόνα υπερέχουν έναντι του BEAV: SN 89.9-93.9% vs 84.8%
Αποτελέσματα διαθέσιμα 24-48 h
νωρίτερα

Επιλογή σύμφωνα με:



ευκολία διάκρισης αποικιών,
αξιοπιστία, κόστος, ανάγκη
επιβεβαιωτικών δοκιμασιών



MALDI-TOF Mass Spectrometry και ανίχνευση αντοχής

Matrix-Assisted Laser Desorption Ionization
(MALDI) Time-Of-Flight (TOF)
Mass Spectrometry



System analysis time to ID result:
95 isolates + 1 QC sample ~5 min






MALDI-TOF MS-based rapid diagnostic methods for AMR

Strategy	Target	Target Resistance	Commercialized (Company) and/or Developing Products	Ref.
Antimicrobial susceptibility testing	AUC ¹ shift spectra (isotope)	All antimicrobials	MBT-ASTRA (Bruker) MBT-Resist (Bruker) DoT-MGA ² (inv. ³)	[22,24,25,70,71]
ID of the AMR clonal group	Clone-specific spectra	Methicillin (MRSA)	ClinProtTools (Bruker)	[72]
ID of the modified antimicrobial drug	Shift peak of the modified drug	Beta-lactams Cephalosporins Carbapenems	MBT STAR [®] -BL (Bruker) MBT STAR [®] -Cepha (Bruker) MBT STAR [®] -Carba (Bruker)	[40,73,74]
ID of the modified antimicrobial target	Shift peak of the modified drug target	Colistin	MALDIxin test (inv.)	[49,75]
Direct detection of the AMR determinant	Protein-specific peak	Carbapenems	Aprodete [™] KPC-CPE (Aprilis)	[55]
Biomarker detection	Protein-specific peak	Carbapenems	ClinProtTools (Bruker)	[64,76]

¹ AUC, area under curve. ² DoT-MGA, Direct-on-Target Microdroplet Growth Assay. ³ inv., under investigation.

MALDI-TOF Mass Spectrometry και ανίχνευση αντοχής

Identification of the KPC producer



beta-lactamases
Porin defects


Direct detection
of AMR determinants

The modified antibiotics
change their molecular
weights.
AMR is identified by
detecting the peak shift
using the MALDI-TOF MS
methodology

Direct detection
of AMR determinants
Identification
of the modified drug

MALDI-TOF MS protein
fingerprints of well-
known AMR
clonal groups

Identification of
the antimicrobial-resistant clonal group
Identification of
the biomarker, non-AMR determinants




MALDI-TOF Mass Spectrometry και ανίχνευση αντοχής

MBT STAR®-Carba IVD Assay




Sensitivity 93.9%
specificity, 100%



Bacteria with active carbapenemase will inactivate the carbapenem antibiotic by hydrolysis of the β -lactam ring, **which is associated with a detectable mass shift.**

In assays with bacteria without active carbapenemase, only peaks corresponding to the intact antibiotic will be present in the mass spectrum.

In assays with bacteria with active carbapenemase, peaks of the intact antibiotic will decrease.



Se: 100%, Sp: 100%

Del Corpo O et al. Clin Microbiol Infect.
2023 Sep 16:S1198-743X(23)00425-1.

Sample	Species	Control ID	Carba
Sample 1	Klebsiella pneumoniae	confirmed	1.09 H
Sample 2	Klebsiella pneumoniae	confirmed	0.98 H
Sample 3	Escherichia coli	confirmed	0.02 NH
neg.control		not performed	-0.03 NH
pos.control		not performed	1.09 H
Δ controls			1.47

H	Hydrolyzed ¹
NH	Non-hydrolyzed ²

MALDI-TOF Mass Spectrometry και ανίχνευση αντοχής

MBT STAR®-Cepha IVD Assay

-  Positive blood culture^s or culture plate
 -  Optional species ID before / in parallel (control ID) with assay*
 -  Bacteria suspension in antibiotic solution
 -  30 min incubation
 -  MALDI target plate preparation
 -  MALDI-TOF MS spectra acquisition and automated data analysis by MBT STAR-BL IVD Module
- 1 hour




Se: 94%

Sp: 97%

Del Corpo O et al. Clin Microbiol Infect.
2023 Sep 16:S1198-743X(23)00425-1.

Hydrolysis of β -lactam ring leads to mass shifts that can easily be detected by MALDI-TOF mass spectrometry




Sample	Species	Control ID	Cepha
Sample 1	Escherichia coli	confirmed	1.28 H
Sample 2	Escherichia coli	confirmed	1.14 H
Sample 3	Klebsiella variicola	confirmed	-0.55 NH
neg.control		not performed	-0.10 NH
pos.control		not performed	1.04 H
Δ controls			1.12

H	Hydrolyzed ¹
NH	Non-hydrolyzed ²

MALDI-TOF Mass Spectrometry και ανίχνευση αντοχής

MBT Subtyping IVD Module

Quick detection of *bla*_{KPC} expressing *K. pneumoniae* and *E. coli*



*A peak in MALDI-TOF mass spectra of KPC-producing *K. pneumoniae*, related to the plasmid carrying *bla*_{KPC} (pKpQIL). This peak at 11,109 m/z is clearly detectable in MALDI-TOF mass spectra.

Prerequisite for the automated detection process is the successful ID.

The MBT HT Subtyping IVD Module then looks for the *bla*_{KPC} related peak in the sample spectrum and, if present, the software will report this sample as a KPC positive one.

KPC detection of *K. pneumoniae* and *E. coli* includes only strains with a *bla*_{KPC} pKpQIL plasmid.

If the gene expression rate is low, there will be no characteristic peak and no KPC subtyping alert.

Cultivation media and conditions might suppress or induce a signal at m/z 11,109 which is not related to a KPC resistance. Therefore, Columbia Agar with 5% sheep blood agar must be used for cultivation of *K. pneumoniae* and *E. coli*.

MALDI-TOF Mass Spectrometry και ανίχνευση αντοχής

MBT Subtyping IVD Module


Bacteroides fragilis *cfiA* subtyping

A cluster of *B. fragilis*, (division II), carries the *cfiA* gene encoding a MBL.

The intrinsic gene is overexpressed when an insertion sequence is integrated upstream of the gene, providing a strong promoter. Then, *cfiA* MBL confers resistance against nearly all beta-lactam drugs to the bacterial host.

Through an extensive screening of *B. fragilis*, the *cfiA*-positive division II strains were distinguished from the *cfiA*-negative strains by spectral shifts between 4000–5500 *m/z*, and the results were satisfactory.

After successful identification, the MBT HT Subtyping IVD Module looks for specific peaks associated with the protein expressed by *B. fragilis* *cfiA* positive strains.



Applications of MALDI-TOF MS for specific AMR detection



Organism	Antibiotic	Year (Reference)
<i>E. coli</i>	Polymyxins	2018 [122]
<i>E. coli</i>	Colistin	2019 [123]
<i>E. coli</i> <i>Klebsiella pneumoniae</i>	Beta-lactams (ESBL-producing isolates)	2019 [124]
<i>Staphylococcus aureus</i> <i>S. intermedius</i> <i>S. pseudintermedius</i>	Novobiocin Polymyxin B Acriflavine	2019 [125]
<i>S. aureus</i>	Methicillin	2019 [126]
<i>Candida auris</i>	Echinocandins	2019 [127]
<i>K. pneumoniae</i> <i>Bacteroides fragilis</i> <i>S. aureus</i>	Carbapenems (carbapenemase-producing isolates) Methicillin	2019 [128]
<i>Enterobacteriaceae</i>	Carbapenems (carbapenemase-producing isolates)	2019 [129]
<i>Enterobacteriaceae</i>	Carbapenem	2019 [130]
<i>Pseudomonas aeruginosa</i>	Beta-lactams (MBL)	2019 [131]
<i>Enterococcus faecium</i>	Vancomycin	2019 [132]
<i>K. pneumoniae</i>	Carbapenems (carbapenemase-producing isolates)	2019 [133]
<i>Acinetobacter baumannii</i>	Colistin	2020 [134]
<i>E. coli</i> <i>K. pneumoniae</i>	Cefotaxime Meropenem, Ciprofloxacin	2020 [135]
<i>S. aureus</i> <i>Enterococcus</i> species <i>E. coli</i> <i>K. pneumoniae</i>	Oxacillin (methicillin) Vancomycin Ceftriaxone Meropenem	2020 [136]
<i>Enterobacterales</i>	Imipenem/Relebactam	2020 [137]
<i>S. aureus</i>	Methicillin	2020 [138]

Kaprou GD et al. Antibiotics (Basel). 2021;10(2):209.



MALDI-TOF MS approaches for AMR detection in *S. aureus*

MALDI-TOF System Used	Principle of the Method	On Culture/ on Clinical Sample	TAT ¹	Brief Advantages/ Disadvantages +/–	No. and Type of <i>S. aureus</i> Analyzed ²	Major Diagnostic Performance ³	Reference
MALDI Biotyper® CA System (Bruker Daltonics)	21 peaks to discriminate MRSA/MSSA	Culture	Few m	+ Use a combination of peaks – Low specificity	181 <i>S. aureus</i>	Sensitivity 87.6% Specificity 71.4%	[64]
MALDI Biotyper® CA System (Bruker Daltonics)	Peptide PSM- <i>mec</i> and δ -toxin peaks for MRSA identification	Culture	Few m	+ Sensitive for SCC <i>mec</i> type II, III and VIII – Not adapted for other SCC <i>mec</i> types	35 MRSA	Sensitivity 40%	[65]
MALDI Biotyper® CA System (Bruker Daltonics)	Peptide PSM- <i>mec</i> peak for MRSA identification	Culture	Few m	+ Good specificity – SCC <i>mec</i> types of isolates unknown	241 MRSA 106 MSSA	Sensitivity 60.2% Specificity 100%	[66]
MALDI Biotyper® CA System (Bruker Daltonics)	m/z 4594 peak for MRSA identification	Culture	Few m	+ Potential MRSA biomarker peak – To be tested in an independent set of data	36 MRSA 31 MSSA	Sensitivity 83.3% Specificity 96.8%	[67]
VITEK MS (BioMérieux)	support vector machine algorithm to identify 38 peaks to discriminate MRSA/MSSA	Culture	Few m	+ Use a combination of peaks – MRSA/MSSA classification model needs optimization	194 MRSA 258 MSSA	Sensitivity 84.0% Specificity 88.0%	[68]
MALDI Biotyper® CA System (Bruker Daltonics)	DOT-MGA	Positive blood cultures	6 h	+ Rapid; tested with cefoxitin – Could be optimized for additional antibiotics	30 MRSA 14 MSSA	Sensitivity 100% Specificity 100%	[69]
VITEK MS (BioMérieux)	DOT-MGA	Culture	5 h	+ Tested with oxacillin – Could be optimized for additional antibiotics	20 MRSA 20 MSSA	Sensitivity 100% Specificity 100%	[70]







Μοριακές μέθοδοι: από την PCR στο mNGS


Gene(s) amplification + amplicons detection

Single / Multiplex Endpoint PCR




Διαδικασία εργώδης
Κίνδυνος επιμολύνσεων

Single / Multiplex Real time PCR





Multiplexing and miniaturization
of the hybridization step

PCR / hybridization σε microarray





Εμπορικά
διαθέσιμες
μέθοδοι με:
αυτοματοποίηση
πολυπλεκτικότητα
σανταρισμένα
πρωτόκολλα
quality control




The basic workflow of molecular-based techniques for AMR detection

Real-time PCR (qPCR)



- measurement of data in real-time
- greater sensitivity
- reduced risk of contamination
- greater amenity to multiplexing
- many systems are partially or even completely automated




Ισοθερμικές μοριακές μέθοδοι και ανίχνευση αντοχής: παράδειγμα

Loop-mediated isothermal amplification technology and a real-time PCR platform

No DNA/RNA extraction necessary

eazyplex[®], Amplex



TAT= 20-30 min

Sample: swab, culture

eazyplex[®] SuperBug CRE

KPC, NDM, OXA-48, OXA-181 and VIM
ESBL genes: CTX-M-1 and CTX-M-9 group

TAT= 15-20 min

Direct screening from **rectal swab and urine**
Confirmation from **positive BC and agar plate**

eazyplex SuperBug Acinetobacter

OXA-51, OXA-23-, -40- and -58-like
and NDM
in 15 minutes

	SuperBug complete A	SuperBug complete B	SuperBug complete C
NDM	X	X	X
VIM	X	X	X
KPC	X	X	X
OXA-48	X	X	X
OXA-23	X	X	
OXA-40	X	X	
OXA-58	X		
OXA-181		X	X
IMP			X

Μοριακές μέθοδοι και ανίχνευση αντοχής: παραδείγματα

Xpert® Carba-R, Cepheid

TAT: 48-minutes



Targeted genes

KPC, NDM, VIM, IMP, OXA-48
(covering OXA-181 & OXA-232)

Δείγμα

Ορθικός
στυλεός
Αποικίες

¹Overall percentage agreement with WGS: 92.8%

²Metanalysis of 17 studies involving 15,972 samples

Pooled Se= 95%, Sp= 99%

In high-risk populations: Se= 99%, Sp= 98%


1. Khoo BY et al. J Clin Microbiol. 2023;61(9):e0031623.

2. Bai Y et al. J Mol Diagn. 2021;23(11):1534-1544.

Check-Direct CPE® assay, (Checks points)



TAT: 2 h
(+DNA extraction)



Targeted genes	Specimen	Technology/ instrumentation
KPC, VIM, NDM, OXA-48	Rectal swab Strains	Two-step ligation-mediated multiplex real-time PCR amplification – micro array hybridization / non-dedicated platform

867 samples from 627 patients (Reference method: culture + PCR)

CPO assay Se= 95.7%, SP= 96.5%, PPV= 60.8% and NPV= 99.8%

Xpert assay

Se= 97.9%, Sp= 99.8%, PPV= 95.8%, and NPV= 99.9%

The CPO assay: useful when handling many specimens, as tests are conducted in batches.

Positive cases showing high Ct values should be confirmed by another assay to rule out false positivity.

Yoo IY et al. Diagn Microbiol Infect Dis. 2022;103(3):115716.

Συμβατική vs. ταχεία διάγνωση



Εμπειρική
αγωγή



Ολικό αίμα



Χρώση Gram



T2Bacteria Panel

Sensitivity: 95%^{1*}
Specificity: 98%^{1*}

E. faecium
S. aureus
K. pneumoniae
P. aeruginosa
E. coli

T2RESISTANCE PANEL

Gram-negative marker
KPC
OXA-48 Group
NDM /VIM/IMP
CTX-M 14/15
AmpC(CMY/DHA)

Gram-positive marker
vanA/B
mecA/C



T2Candida Panel

Sensitivity: 91%³
Specificity: 99%³

Candida albicans
Candida tropicalis
Candida krusei
Candida glabrata
Candida parapsilosis

Ροή εργασίας συμβατικών μεθόδων (~48-168 h)

Short subculture → ID /AST



Ροή εργασίας
μοριακών
μεθόδων

Διάγνωση
Στοχευμένη
αγωγή





Θετική

BC

Direct MALDI-TOF MS identification from positive blood cultures

ADVANTAGES

Very rapid identification
Easy
Low additional cost

LIMITATIONS

Identification difficulties with some species
Confirmation might be needed
Mixed cultures

MALDI-TOF MS identification from positive blood cultures using short sub-culture on solid medium

ADVANTAGES

Rapid identification
Very easy
Easy to integrate in lab workflow
Confirmation usually not needed

LIMITATIONS

Slow-growing organisms
Mixed cultures

Molecular methods for identification from positive blood cultures

ADVANTAGES

Rapid identification
Easy

LIMITATIONS

High cost
Confirmation needed
Limited pathogen panels

Reduced diagnostic accuracy for polymicrobial cultures

Sensitivity is influenced by:

- the concentration of pathogen-derived nucleic acids
- the presence of PCR inhibitors
- the sample volume

DNA-based identification from whole blood

ADVANTAGES

Time advantage compared to culture
Advantage in pathogen detection under antibiotics

LIMITATIONS

Limited number of pathogens in panel assays
Relatively low sensitivity (only 1-5 ml blood sampled)
Resistance detection ≠ susceptibility testing
Clinical relevance of DNAemia ?
High workload
Highly experienced staff required
High cost
Analytical cut-offs set by some manufacturers
Does not substitute culture-based testing

90% των αιτίων μικροβιαμίας

Ολικό αίμα



Prevalence of Antibiotic-Resistant Pathogens in Culture-Proven Sepsis and Outcomes Associated With Inadequate and Broad-Spectrum Empiric Antibiotic Use

Resistant GP isolated in 13.6% of patients

Resistant GN in 13.2%

Both undertreatment (failure to cover organisms)

and overtreatment

(resistant organisms targeted but not isolated)

were associated with higher mortality after detailed risk adjustment.

N= 17 430 adults with culture-positive sepsis

admitted to 104 US hospitals

67% received empiric broad-spectrum antibiotics

Is there risk associated with the receipt of empiric broad-spectrum antibiotics?

Table 2.

Outcomes Associated With Inadequate and Unnecessarily Broad Empiric Antibiotic Therapy^a

Outcome	Inadequate vs adequate empiric therapy						Unnecessarily broad vs not unnecessarily broad empiric therapy ^b					
	No./total No. (%)		Unadjusted OR (95% CI)	P value	Adjusted OR (95% CI)	P value	No./total No. (%)		Unadjusted OR (95% CI)	P value	Adjusted OR (95% CI)	P value
	Inadequate	Adequate empiric therapy					Unnecessarily broad	Not unnecessarily broad				
In-hospital death	488/2785 (17.5)	2011/12 388 (16.3)	1.10 (0.98-1.22)	.09	1.19 (1.03-1.37)	.02	1575/8405 (18.7)	436/3993 (10.9)	1.88 (1.68-2.11)	<.001	1.22 (1.06-1.40)	.007
Hospital-onset acute kidney injury	486/2785 (17.5)	2196/12 398 (17.7)	0.98 (0.88-1.09)	.74	1.02 (0.90-1.16)	.72	1641/8405 (19.5)	555/3993 (13.9)	1.50 (1.35-1.67)	<.001	1.12 (1.00-1.26)	.05
<i>Clostridioides difficile</i>	207/2785 (7.4)	498/12 398 (4.0)	1.92 (1.63-2.27)	<.001	1.19 (0.98-1.45)	.09	367/8405 (4.4)	131/3993 (3.3)	1.34 (1.10-1.65)	.004	1.26 (1.01-1.57)	.04

These findings underscore the need for better tests to rapidly identify patients with resistant pathogens

FilmArray BCID2: 43 Targets, TAT 1 h

GRAM-NEGATIVE BACTERIA

Acinetobacter baumannii complex

Bacteroides fragilis

Enterobacteriales

Enterobacter cloacae complex

Escherichia coli

Klebsiella aerogenes

Klebsiella oxytoca

Klebsiella pneumoniae group

Proteus spp.

Salmonella spp.

Serratia marcescens

Haemophilus influenzae

Neisseria meningitidis

Pseudomonas aeruginosa

Stenotrophomonas maltophilia

GRAM-POSITIVE BACTERIA

•*Enterococcus faecalis*

•*Enterococcus faecium*

•*Listeria monocytogenes*

•*Staphylococcus* spp.

Staphylococcus aureus

Staphylococcus epidermidis

Staphylococcus lugdunensis

•*Streptococcus* spp.

Streptococcus agalactiae

Streptococcus pneumoniae

Streptococcus pyogenes

YEAST:

•*Candida albicans*

•*Candida auris*

•*Candida glabrata*

•*Candida krusei*

•*Candida parapsilosis*

•*Candida tropicalis*

•*Cryptococcus*

(*C. neoformans*/*C. gattii*)

AMR GENES

Carbapenemases

IMP

KPC

OXA-48-like

NDM

VIM

Colistin Resistance: *mcr-1*

ESBL: CTX-M

Methicillin Resistance


mecA/C

mecA/C and MREJ


(MRSA)

Vancomycin Resistance:

vanA/B



5 Inject Sample into Pouch



6 Insert Pouch into FilmArray and Start Run

Nested-PCR and melting curve analysis dedicated platform

ePlex® BCID Panels



*Se 97-99%

Sp 96-99%

*Huang TD et al.

J Clin Microbiol. 2019;57(2).

Gram-positive Organisms

21 παθογόνα
6 γονίδια αντοχής
Gram-Negative Organisms

20 παθογόνα
4 γονίδια αντοχής

15
παθογόνα

Fungal Organisms

VERIGENE® Bloodstream Infection Testing Panels

Gram-Positive Blood Culture Test (BC-GP)

Species

Staphylococcus aureus
Staphylococcus epidermidis
Staphylococcus lugdunensis
Streptococcus agalactiae
Streptococcus pneumoniae
Streptococcus pyogenes
Enterococcus faecalis
Enterococcus faecium

Group

Streptococcus anginosus

Genus

Staphylococcus spp.
Streptococcus spp.
Micrococcus spp.*
Listeria spp.

Resistance

mecA (methicillin)
vanA (vancomycin)
vanB (vancomycin)

Gram-Negative Blood Culture Test (BC-GN)

Species

*Escherichia coli**
Klebsiella pneumoniae
Klebsiella oxytoca
Pseudomonas aeruginosa
*Serratia marcescens***

Genus

Acinetobacter spp.
Citrobacter spp.
Enterobacter spp.
Proteus spp.

Resistance

CTX-M (ESBL)
IMP (carbapenemase)
KPC (carbapenemase)
NDM (carbapenemase)
OXA (carbapenemase)
VIM (carbapenemase)

Se: 98.1% GN, 97.3% GP, 99.2% *Candida*

Sp: 99.9% GN, 99.8% GP, 99.9% *Candida*

Se/Sp για AMR γονίδια: 98.4-100/98.3-100%

Salimnia H et al. J Clin Microbiol 2016.

Micro-array hybridization dedicated platform



TAT 2-2.5 h



Συνδρομική διάγνωση και ανίχνευση γονιδίων αντοχής: Blood Culture Panels (positive BC bottle)

Parameter	FilmArray BCID	Verigene	
		Gram-positive blood culture	Gram-negative blood culture
Total no. of targets	27	15	14
Ability to detect presence of resistance gene			
<i>mecA</i>	✓	Se 93.9%	✓
<i>vanA</i>	✓	Se: 90.9	✓
<i>vanB</i>	✓		
<i>bla</i> _{KPC}	✓		
<i>bla</i> _{NDM}	✓		
<i>bla</i> _{OXA}	✓		
<i>bla</i> _{VIM}	✓		
<i>bla</i> _{IMP}	✓		
<i>bla</i> _{CTX-M}	✓		
<i>mcr-1</i>	✓		
Time to result (h)	~1	~2.5	~2

Sp > 98 %

Se 93.9%


Se: 90.9

NPA >99.9%

Se: 94.9%

PPA 100%
PPA 96.2%
PPA 94.3%
PPA 100%
PPA 100%
PPA 98.9%

ePlex BCID, GenMark



Multiplexed NAAT and a electrochemical detection technology with a dedicated platform

Peri AM et al. BMC Infect Dis. 2022;22(1):794.




Ramanan P et al. Clin Microbiol Rev. 2017;31(1) (με επικαιροποίηση).

Salimnia H et al. J Clin Microbiol. 2016;54(3):687-98.

Ledeboer NA et al. J Clin Microbiol. 2015;53(8):2460-72.

Prospective RCT, 617 patients with BSI

Gram result reported






ID Organism identification
 AST Phenotypic antimicrobial susceptibility report
 ▽ D De-escalation
 △ E Escalation

2022;75(12):2066-2075.

Impact of a Rapid Molecular Test for *Klebsiella pneumoniae* Carbapenemase and Ceftazidime-Avibactam Use on Outcomes After Bacteremia Caused by Carbapenem-Resistant Enterobacterales

Michael J. Satlin,^{1,2} Liang Chen,^{3,4} Angela Gomez-Simmonds,⁵ Jamie Marino,² Gregory Weston,⁶ Tanaya Bhowmick,⁷ Susan K. Seo,⁸ Steven J. Sperber,^{9,10} Angela C. Kim,¹¹ Brandon Eilertson,¹² Sierra Derti,¹ Stephen G. Jenkins,^{1,2} Michael H. Levi,¹³ Melvin P. Weinstein,^{7,14} Yi-Wei Tang,¹⁵ Tao Hong,¹⁶ Stefan Juretschko,¹⁷ Katherine L. Hoffman,¹⁸ Thomas J. Walsh,¹ Lars F. Westblade,^{1,2} Anne-Catrin Uhlemann,⁵ and Barry N. Kreiswirth^{3,4}

N= 137 patients with CRE bacteremia**[77% CP-CRE (65% KPC)]****2016-2018, 8 New York and New Jersey medical centers****3/8 study hospitals used the FilmArray BCID**

30-day mortality
24% vs 47%;
 $P=.007$

30-day
mortality 10% with
caz-avi monotherapy
vs. 31% with
polymyxin
monotherapy
($P=.08$)

Impact of an Antimicrobial Stewardship Program Intervention Associated with the Rapid Identification of Microorganisms by MALDI-TOF and Detection of Resistance Genes in ICU Patients with Gram-Negative Bacteremia

216 episodes of GN bacteremia, ICU

MALDI-TOF και ανίχνευση των γονιδίων AMR απευθείας από τη φιάλη της θετικής αιμοκ/ας

Table 4. Antimicrobial consumption in DOT per 1000 days present and collection (21 days).

Antimicrobial Consumption and Cost	Pre-Intervention	Intervention	p-Value
Consumption (DOT/1000 days present, median, IQR)			
All antimicrobials	N = 114 1.381 (1.103–2.251)	N = 102 1.262 (1.063–1.662)	0.032 ^a
Antimicrobial for gram-negative bacteria	N = 114 1.281 (1.004–1.775)	N = 102 1.172 (1.006–1.427)	0.067 ^a
Antimicrobial for gram-positive bacteria	N = 52 475 (238–761)	N = 43 270 (139–467)	0.004 ^a
Carbapenems	N = 76 836 (504–1056)	N = 72 543 (301–991)	0.040 ^a
Colistin/Polymyxin B	N = 45 722 (390–932)	N = 40 881 (462–1001)	0.299 ^a
Ceftazidime/avibactam	N = 3 931 (725–1022)	N = 8 911 (823–940)	0.865 ^a

Καμιά επίδραση στη θνητότητα 30-ημερών
Σημαντικά μικρότερη η συνολική διάρκεια νοσηλείας (44 vs. 39 ημέρες, $p = 0.005$) και ο χρόνος νοσηλείας στη ΜΕΘ (17 vs. 13 ημέρες, $p = 0.033$)



EXPERT REVIEW OF MEDICAL DEVICES
2021, VOL. 18, NO. 5, 473-482
<https://doi.org/10.1080/17434440.2021.1919508>

META-ANALYSIS

OPEN ACCESS Check for updates

Antimicrobial and resource utilization with T2 magnetic resonance for rapid diagnosis of bloodstream infections: systematic review with meta-analysis of controlled studies



61.4% των αιτίων
μικροβιαμίας



T2BACTERIA PANEL

Sensitivity: 90%⁷
Specificity: 98%⁷

E. faecium
S. aureus
K. pneumoniae
P. aeruginosa
E. coli



T2RESISTANCE PANEL

Gram-negative marker
KPC
OXA-48 Group
NDM /VIM/IMP
CTX-M 14/15
AmpC(CMY/DHA)

Gram-positive marker
vanA/B
mecA/C




T2CANDIDA PANEL

Sensitivity: 91%⁵
Specificity: 99%⁵

Candida albicans
Candida tropicalis
Candida krusei
Candida glabrata
Candida parapsilosis


1-3 CFU/mL LoD

Hours to species identification with T2MR vs. BC




Pooled mean difference = -77 hours (p < 0.001)

Hours to targeted therapy among T2MR positive cases vs. BC




Pooled mean difference = -42 hours (p < 0.001)

Hours to empirical therapy de-escalation among T2MR negative cases vs. BC




Pooled mean difference = -7 hours (p = 0.02).

Mortality with T2MR vs. BC




28.9% vs. 29.9%, RR = 1.02, p = 0.86

ICU stay with T2MR vs. BC. Pooled mean difference = -5.0 days (p =0.03)



Hospital stay with T2MR vs. BC. Pooled mean difference = -4.8 days (p =0.04)



Considering hospital costs for sepsis patients range from 2,000 USD to 5,000 USD per day depending on sepsis severity, implementation of T2MR could theoretically reduce hospital costs by 10,000 USD to 25,000 USD per patient tested

Μοριακή ανίχνευση αντοχής: προβληματισμοί

Ασυμφωνία αποτελεσμάτων μοριακών / συμβατικών μεθόδων

Ποια είναι η μέθοδος αναφοράς / gold standard;

Υψηλά ποσοστά **culture-negative / molecular testing-positive** (5% - 55%)

- Αυξανόμενη ευαισθησία σύγχρονων μοριακών μεθόδων
Πιθανή έλλειψη ευαισθησίας καλλιεργητικών μεθόδων (low inoculum)
- Παρουσία νεκρών ή μη καλλιεργήσιμων βακτηρίων
- Μη έκφραση του γονιδίου αντοχής ή έκφραση σε χαμηλά επίπεδα
(χαμηλού επιπέδου αντοχή μη ανιχνεύσιμη με τις συμβατικές μεθόδους)
- Παρουσία άλλων βακτηρίων στο δείγμα που δεν ανιχνεύονται
και φέρουν τα γονίδια αντοχής



Ποιες είναι οι κλινικές επιπτώσεις


FP μοριακές, ME (μη απαραίτητη ευρεία κάλυψη ασθενούς)

FN καλλιεργητικές, VME ακατάλληλη αγωγή

Περιορισμένες επιπτώσεις στην περίπτωση μη έκφρασης του γονιδίου

Μοριακή ανίχνευση αντοχής: προβληματισμοί

Detection of particular resistance mechanisms vs. phenotypic AST



- 
- Πόσα και ποια γονίδια περιλαμβάνονται στα panel
 - Επιλογή γονιδίων σε σχέση με την τοπική επιδημιολογία τη δεδομένη χρονική στιγμή
 - Μη αναγνώριση κάποιων variants ή αναγνώριση variants με διαφορετικό υδρολυτικό προφίλ
 - Μη ανίχνευση νέων αναδυόμενων γονιδίων αντοχής
 - **Σύνθετοι μηχανισμοί αντοχής**
 - **Πολυμικροβιακά δείγματα και γονίδια αντοχής**

Είναι ασφαλής η αποκλιμάκωση;



Δεν υπάρχει σαφής κατεύθυνση για τις περιπτώσεις που δεν ανιχνεύονται καθοριστές αντοχής

An Antibiotic Stewardship Program Blueprint for Optimizing Verigene BC-GN within an Institution: a Tale of Two Cities

1046 GN bloodstream isolates

TABLE 2 Verigene BC-GN analytic results by bug-drug combination^a

Institution and target organism	n	No. (%) of isolates	No. of isolates identified by	P (%)	PPV (%)	NPV (%)
Detroit Medical Center						
<i>E. coli</i>	31	9	97	98		
<i>K. pneumoniae</i>	14	9	97	94		
	14	00	100	99		
<i>K. oxytoca</i>	21	00	100	95		
<i>Proteus</i> spp.	5	00	100	94		
<i>Enterobacter</i> spp.	6	00	100	97		
<i>Citrobacter</i> spp.	10			100		
<i>Acinetobacter</i> spp.	39	3	80	93		
<i>Pseudomonas</i> spp.	5			88		



AMERICAN SOCIETY FOR MICROBIOLOGY

Journal of Clinical Microbiology®

2020;58(4):e02098-19.

LETTER TO THE EDITOR

Check for updates

The Clinical Impact of a Negative Molecular β -Lactamase Gene Test for *Enterobacteriaceae*: Let's Not Let Perfect Be the Enemy of Really Good

Jason M. Pogue,^a Emily L. Heil^b

Με εξαίρεση την *P. aeruginosa*, η απουσία καθοριστών αντοχής μπορούσε να προβλέψει την ευαισθησία στα αντιβιοτικά στόχους



ID and phenotypic AST using the Accelerate Pheno System (RAPID)

RCT

Patients (N= 448)

with positive BC with GNB

RANDOMIZED

SOC testing with AS review
or RAPID with AS





Patients with MDRO, escalation

Median RAPID, 18.4 vs SOC, 61.7 h
(difference, 43.3 h; $P = .01$)

Patients with S organisms, de-escalation

median RAPID, 29.4 vs SOC, 42.4 h
(difference, 13 h, $P = .02$)

Banerjee R et al. Clin Infect Dis. 2021;73(1):e39-e46.



Real-World Impact of the Accelerate PhenoTest BC Kit on Patients With Bloodstream Infections in the Improving Outcomes and Antimicrobial Stewardship Study: A Quasiexperimental Multicenter Study

Amira A. Bhalodi,^{1a} Shawn H. MacVane,^{1a} Bradley Ford,² Dilek Ince,³ Patrick M. Kinn,⁴ Kelly M. Percival,⁴ Derek N. Bremmer,⁵ Dustin R. Carr,⁵ Thomas L. Walsh,⁶ Micah M. Bhatti,⁷ Samuel A. Shelburne,⁸ Romney M. Humphries,^{1b} Kaleb Wolfe,⁹ Eric R. Rosenbaum,¹⁰ Ryan K. Dare,⁹ Johann Kolev,¹¹ Meghan Madhusudhan,¹² Michael A. Ben-Aderet,¹¹ and Margie A. Morgan¹³

N= 854 ασθενείς με BSIs
5 κέντρα στις ΗΠΑ
Απρίλιος 2017 - Νοέμβριος 2019





Median TTOT was 17.2 hours shorter in the post-AXDX arm (23.7 h) vs. pre-AXDX arm (40.9 h; $P < .0001$).


Διάμεσος χρόνος

Για την πρώτη τροποποίηση: 24.2 vs 13.9 h; $P < .0001$

Για την πρώτη αποκλιμάκωση: 36.0 vs 27.2 h; $P = .0004$



Rapid AST



RAPID AST SYSTEM 5.5 hours



TABLE 1 Overall performance of the Reveal AST system

Parameter or detail	Performance of Reveal AST against:	
	Sensititre	Vitek 2
Parameter, % (no. positive/total no.)		
EA ^c	98.0 (2,129/2,173)	97.0 (1,482/1,528)
CA ^c	96.3 (2,174/2,258)	96.2 (1,554/1,615)
mE	3.5 (78/2,258)	3.3 (54/1,615)
ME	0.3 (5/1,889)	0.3 (4/1,342)
VME	1.3 (4/313)	1.3 (3/232)


Wide antimicrobial coverage for Gram-negative bloodstream infections, with 176 bug drug combinations.

ANTIMICROBIAL DRUG	MIC CALLING RANGE (µg/ml)	A. baumannii	C. freundii	C. koseri	E. cloacae	E. coli	K. aerogenes	K. oxytoca	K. pneumoniae	P. mirabilis	P. aeruginosa	
Amikacin	4 - 16	✓	✓	✓	✓	✓	✓	✓	✓	✓	10	
Amoxicillin / Clavulanate	4/2 - 16/2			✓	✓		✓	✓			5	
Ampicillin	4 - 8				✓				✓		2	
Aztreonam	1 - 16		✓	✓	✓	✓	✓	✓		✓	8	
Cefepime	0.125 - 64		✓	✓	✓	✓	✓	✓	✓	✓	8	
Cefotaxime	0.125 - 4		✓	✓	✓	✓	✓	✓	✓		8	
Cefoxitin	8 - 16			✓	✓		✓	✓	✓		5	
Ceftazidime	0.125 - 64		✓	✓	✓	✓	✓	✓	✓	✓	9	
Ceftazidime / Avibactam	0.25/4, 1/4 - 16/4		✓	✓		✓	✓	✓	✓	✓	9	
Ceftazidime / Clavulanate	0.25/4 - 8/4				✓		✓	✓			3	
Certolozane / Tazobactam	1/4 - 4/4				✓	✓	✓	✓	✓	✓	6	
Ciprofloxacin	0.06, 0.25 - 1	✓	✓	✓	✓	✓	✓	✓	✓	✓	10	
Ertapenem	0.125 - 1		✓	✓	✓	✓	✓	✓	✓		8	
Gentamicin	2 - 4	✓	✓	✓	✓	✓	✓	✓	✓		9	
Imipenem	1 - 8	✓	✓	✓	✓	✓	✓	✓	✓	✓	9	
Levofloxacin	0.25 - 1	✓	✓	✓	✓	✓	✓	✓	✓	✓	10	
Meropenem	0.125 - 8	✓	✓	✓	✓	✓	✓	✓	✓	✓	10	
Meropenem / Vaborbactam	2/8 - 8/8		✓	✓	✓	✓	✓	✓	✓	✓	9	
Piperacillin	8 - 16		✓		✓	✓	✓	✓	✓	✓	8	
Piperacillin / Tazobactam	4/4 - 16/4		✓	✓	✓	✓	✓	✓	✓	✓	9	
Tigecycline	0.5 - 1			✓		✓					2	
Tobramycin	2 - 4	✓	✓	✓	✓	✓	✓	✓	✓	✓	10	
Trimethoprim / Sulfamethoxazole	2/38 - 4/76	✓	✓	✓	✓	✓	✓	✓	✓		9	
TOTAL		8	17	19	18	23	17	21	21	18	14	176

System reports results with S/I/R interpretation according to CA-SFM 2019 and EUCAST 2021 Breakpoints

Tibbetts R et al. J Clin Microbiol. 2022;60(6):e0009822.

Mortality outcomes with molecular rapid diagnostic testing (mRDT) versus conventional testing in bloodstream infection



mRDT και ASP

31 μελέτες, 5920 BSI pts

Mortality risk decreased significantly with mRDT in the presence of ASP

Number needed to treat: 20 patients with BSI to prevent 1 death within 30 days

Time to effective therapy decreased by a weighted mean difference of –5.03 hours

Length of stay decreased by –2.48 days





The Cost-Effectiveness of Rapid Diagnostic Testing for the Diagnosis of Bloodstream Infections with or without Antimicrobial Stewardship

Elina Eleftheria Pliakos,^a Nikolaos Andreatos,^a Fadi Shehadeh,^a Panayiotis D. Ziakas,^a Eleftherios Mylonakis^a




TABLE 3 Base-case analysis results for competing strategies

Strategy	Base-case estimate				
	Cost (\$)	Effect		ICER	
		Probability of survival	QALY value	\$/death averted	\$/QALY gained
mRDT (with and without ASP)	36,301.50	0.89	11.9883	-490,763	-36,434
mRDT with ASP (A)	31,274.24	0.89	11.9883	-616,445	-45,764
mRDT without ASP (B)	57,220.14	0.90	12.1230	25,762	1,913
PCR with ASP (A1)	47,917.58	0.88	11.8536	-267,148	-19,833
MALDI-TOF analysis with ASP (A2)	28,394.21	0.92	12.3924	-393,397	-29,205
PNA-FISH with ASP (A3)	53,226.09	0.85	11.4495	0	0
BC-GP with ASP (A4)	24,904.91	0.84	11.3148	Dominated	Dominated
BC-GN with ASP (A5)	33,691.47	0.92	12.3924	-317,722	-23,587
PCR without ASP (B1)	47,430.21	0.88	11.8536	-283,394	-21,039
PNA-FISH without ASP (B2)	58,284.16	0.92	12.3924	33,602	2,495
Conventional laboratory methods with ASP (C)	41,723.98	0.84	11.3148	Dominated	Dominated
Conventional laboratory methods without ASP (D)	55,932.02	0.85	11.4495	Baseline	Baseline



Ταχεία μοριακή διάγνωση και diagnostic stewardship

Right test



Right patient



Right time



TABLE 1 Key diagnostic stewardship considerations for implementation of rapid infectious disease diagnostics

Goal	Key question	Key considerations and potential strategies
Right test	Is the test appropriate for the clinical setting?	Sensitivity and specificity Predictive values Testing volumes Diagnostic yield Laboratory feasibility Cost Clinical impact
Right patient	Will the clinical care of the patient be affected by the test result?	Laboratory test utilization committee Automatic laboratory reflex CPOE decision support Appropriate use criteria Indication selection Prior authorization Benchmarking Specimen rejection
Right time	Will the result be available in time to optimally affect care?	Time to specimen receipt Centralized vs point-of-care testing On-demand vs batched testing Specimen preparation time Run time Result reporting time

Λειτουργεία 24/7?

Τρόπος παρουσίασης αποτελεσμάτων

Μοριακή ανίχνευση: αναφορά αποτελέσματος

Reporting

“detected” or “not detected” for a specific molecular target

+

interpretative comments with therapeutic guidance

Simner PJ et al. Clin Infect Dis. 2023;76(9):1550-1558.

Organisms Detected: Enterobacterales Klebsiella pneumoniae group		Controls: Passed	
Applicable Antimicrobial Resistance Genes Detected: CTX-M NDM			
Note: Antimicrobial resistance can occur via multiple mechanisms. A Not Detected result for antimicrobial resistance gene(s) does not indicate antimicrobial susceptibility. Subculturing is required for species identification and susceptibility testing of isolates.			
Result Summary			
Gram Negative Bacteria			
✓	Not Detected	Acinetobacter calcoaceticus-baumannii complex	
	Not Detected	Bacteroides fragilis	
	Detected	Enterobacterales	
	Not Detected	Enterobacter cloacae complex	
	Not Detected	Escherichia coli	
	Not Detected	Klebsiella aerogenes	
	Not Detected	Klebsiella oxytoca	
✓	Detected	Klebsiella pneumoniae group	
	Not Detected	Proteus spp.	
	Not Detected	Salmonella spp.	
	Not Detected	Serratia marcescens	
	Not Detected	Haemophilus influenzae	
	Not Detected	Neisseria meningitidis	
	Not Detected	Pseudomonas aeruginosa	
	Not Detected	Stenotrophomonas maltophilia	
Yeast			
	Not Detected	Candida albicans	
	Not Detected	Candida auris	
	Not Detected	Candida glabrata	
	Not Detected	Candida krusei	
	Not Detected	Candida parapsilosis	
	Not Detected	Candida tropicalis	
	Not Detected	Cryptococcus neoformans/gattii	

Organisms Detected: Enterobacterales Klebsiella pneumoniae group		Controls: Passed	
Applicable Antimicrobial Resistance Genes Detected: CTX-M NDM			
Note: Antimicrobial resistance can occur via multiple mechanisms. A Not Detected result for antimicrobial resistance gene(s) does not indicate antimicrobial susceptibility. Subculturing is required for species identification and susceptibility testing of isolates.			
Result Summary			
Antimicrobial Resistance Genes			
✓	Detected	CTX-M	
	Not Detected	IMP	
	Not Detected	KPC	
	Not Detected	mcr-1	
⊖	N/A	mecA/C	
⊖	N/A	mecA/C and MREJ (MRSA)	
✓	Detected	NDM	
⊖	Not Detected	OXA-48-like	
	N/A	vanA/B	
	Not Detected	VIM	
Gram Positive Bacteria			
	Not Detected	Enterococcus faecalis	
	Not Detected	Enterococcus faecium	
	Not Detected	Listeria monocytogenes	
	Not Detected	Staphylococcus spp.	
	Not Detected	Staphylococcus aureus	
	Not Detected	Staphylococcus epidermidis	
	Not Detected	Staphylococcus lugdunensis	
	Not Detected	Streptococcus spp.	
	Not Detected	Streptococcus agalactiae (Group B)	
	Not Detected	Streptococcus pneumoniae	
	Not Detected	Streptococcus pyogenes (Group A)	

Μοριακή ανίχνευση αντοχής: σχόλια για την αναφορά των αποτελεσμάτων

Table 3. Reporting comments for molecular blood culture diagnostics.

Adapted from Banerjee et al.²⁹


<i>Staphylococcus aureus</i> , <i>mecA</i> detected	Probable methicillin-resistant <i>Staphylococcus aureus</i> (MRSA); further testing in progress. MRSA is predictably resistant to beta-lactam antibiotics (except ceftaroline). Patient requires contact precautions if hospitalized.	<i>Enterobacter cloacae</i> complex	This organism may contain an inducible β -lactamase. Penicillin or second- or third-generation cephalosporin monotherapy may result in emergence of high-level resistance.
<i>Staphylococcus aureus</i> , <i>mecA</i> not detected	Methicillin (oxacillin)-susceptible <i>Staphylococcus aureus</i> . Preferred therapy is an anti-staphylococcal β -lactam antibiotic, unless clinically contraindicated.	<i>Escherichia coli</i> , <i>bla</i> _{CTX-M} detected	Extended-spectrum β -lactamase-producer. Carbapenems are drugs of choice for ESBL-producers.
<i>Staphylococcus</i> , coagulase-negative, <i>mecA</i> detected	Methicillin (oxacillin)-resistant coagulase-negative <i>Staphylococcus</i> . Possible blood culture contaminant (unless isolated from more than one blood culture draw or clinical case suggests pathogenicity). No antibiotic treatment is indicated for blood culture contaminants.	<i>E. coli</i> , <i>bla</i> _{KPC} detected	Carbapenemase producer. Patient requires contact precautions if hospitalized. This organism is resistant to carbapenems and other β -lactam antibiotics. Consult infectious diseases.

The Genotype-to-Phenotype Dilemma: How Should Laboratories Approach Discordant Susceptibility Results?


Rebecca Yee,^a Jennifer Dien Bard,^{ab} Patricia J. Simner^c

difficult to associate the AMR marker with a particular organism
when multiple Gram-negative organisms are identified

Potential steps and reporting structure for specific discrepant examples found in GN




WGS: τα βήματα της διαδικασίας



Metagenomic NGS (agnostic NGS): το επόμενο βήμα

hypothesis-free,
culture-independent,
pathogen detection

evaluate the immunologic
response associated
with the presence and type of infection
e.g. biomarkers
provide insight into the pathogenesis
of the microorganisms detected



Ερμηνεία αποτελεσμάτων

WGS:
“one test fits all”
approach



Tweet versus...

Leo Tolstoy, War and Peace


Napoleon invades Russia.

Russian aristocratic families sent into a tizz


War ensues. French retreat.

Russians celebrate. Lots of them marry.

Tome



Reed T. ECCMID, 2016.



Woodford N. ECCMID, 2016.



Rebelo AR et al. 2022;13:804627.


One Day in Denmark: Comparison of Phenotypic and Genotypic Antimicrobial Susceptibility Testing in Bacterial Isolates From Clinical Settings

N= 500 isolates randomly selected

concordance of 91.7%

Percentage and number of genotype–phenotype concordances and discordances observed for all bacteria analyzed in this study, when comparing MIC and WGS results.

cases of phenotypically resistant isolates without any known genetic resistance mechanism: 2.1% of all combinations analyzed



One size does not fit all



stewardship

Noun: the activity or job of
protecting
and being **responsible**
for something



Our time with antibiotics
is running out.

CHANGE CAN'T WAIT



THE FUTURE OF
ANTIBIOTICS
DEPENDS ON ALL OF US



World AMR
Awareness
Week 2025
18-24 Nov

World AMR Awareness Week

World AMR Awareness Week

18 - 24 November 2025



WAAW 2025 theme

**Act Now:
Protect Our Present,
Secure Our Future**