

Biotechnology course

Polymerase Chain Reaction (PCR)

Dr. Stella Georgiou

sgeorgiou@bioacademy.gr

Central Dogma of Molecular Biology

by
FRANCIS CRICK



MRC Laboratory of Molecular Biology,
Hills Road,
Cambridge CB2 2QH

The central dogma of molecular biology deals with the detailed residue-by-residue transfer of sequential information. It states that such information cannot be transferred from protein to either protein or nucleic acid.


"The central dogma, enunciated by Crick in 1958 and the keystone of molecular biology ever since, is likely to prove a considerable over-simplification."

This quotation is taken from the beginning of an unsigned article¹ headed "Central dogma reversed", recounting the very important work of Dr Howard Temin² and others³ showing that an RNA tumour virus can use viral RNA as a template for DNA synthesis. This is not the first time that the idea of the central dogma has been revised.

analogous to the symbols for
The principle of formulation from one point of view. This could be Fig. 1 (which is not surprising since



multi-OMICS approach



multi-OMICS approach

A genome is an org
Each genome contain
genome—more than


A transcriptome is
The term "transcripto
type. In contrast with

A proteome is the
Proteomics involves t

The metabolome is
during metabolism
The metabolome pro

The epigenome co
gene expression.
The epigenome is a n

Personal "Omics" Profiling (POP)




pcr

pcr

pcr

pcr

pcr



End-point
PCR

RT-PCR

Digital
PCR

Microfluidic
PCR


qRT-PCR

Questions

What is the goal of the experiment?


- **What biological question** is to be answered?
- What is the total number of **genes** to be analyzed?
- What **control samples and genes** will be used to measure the changes in expression levels?
- Are there any **limitations** to the amount of target material available?
- What is the **sensitivity** required to obtain the data necessary to answer?

When and why to check for RNA expression?






- Gene expression changes between different cell types, developmental changes, treatments, diseases/infections/cancer, stress, alternative splicing, etc...

Reverse transcription:




How do we make copies of DNA?

Polymerase Chain Reaction (PCR)



Polymerase Chain Reaction (PCR)

PCR relies on the complementary base pairing of the DNA



PCR uses DNA hybridization

- DNA hybridization involves the formation of a double-stranded nucleic acid
- DNA occurs as a double strand where each *A* is paired with a *T* and vice versa, and each *C* is paired with a *G* and vice versa

The reverse complement of a DNA sequence is formed by reversing the letters, interchanging *A* and *T* and interchanging *C* and *G*

DNA sample >

ATG CCC GGG TAA

Complement : **A** becomes **T**


Reversing : 1st position becomes last

Reverse complement >


T

PCR steps


1. Denaturation




2. Annealing



3. Elongation



Repeat



Fold expression = 2^n
 n = number of cycles





Figure 6. The exponential amplification of DNA in PCR

PCR experiment

What really happens..

How everything starts..



PCR experiment

what usually happens in the lab..



PCR experiment

Successful PCR: everything is optimized

- **Template DNA:** genomic DNA (gDNA), complementary DNA (cDNA), and plasmid DNA, e.g. 5-50 ng of gDNA
- **Buffer:** suitable chemical environment for activity of DNA polymerase. Buffer pH 8.0-9.5
- **DNA polymerase:** in a typical 50 μL reaction, 1-2 units of DNA polymerase
- **Mg²⁺, other cofactors:** Mg²⁺ functions as a cofactor for activity of DNA polymerases
- **Deoxynucleoside triphosphates (dNTPs):** dATP, dCTP, dGTP, and dTTP are in excess. Final concentration of each dNTP is generally 0.2 mM
- **Primers...?**


1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16



← 587 bp
← 500 bp
← 376 bp
← 252 bp
← 211 bp

PCR FW and RV primers

- The **forward (FW)** primer will anneal with 3'-5' DNA strand
- The **reverse (RV)** primer will anneal with the 5'-3' DNA strand




How to choose FW & RV primers?

primers should generally have the following properties:

- Length of 18-24 bases
- 40-60% G/C content
- Start and end with 1-2 G/C pairs
- Melting temperature (T_m) of 50-60°C
- Primer pairs should have a T_m within 5°C of each other
- Primer pairs should not have complementary regions

Note: If you will be including a restriction site at the 5' end of your primer, note that a 3-6 base pair "clamp" should be added upstream in order for the enzyme to cleave efficiently (e.g. GCGGCG-restriction site-your sequence).



Important for the analysis

*Gene-specific!
Or transcript-specific*




PCR products visualization

- In agarose gels the samples are separated by their size
- Electrophoresis is the motion of charged particles relative to a fluid caused by an electric field
- As particle velocities depend on particle **net charge and size**, it can be used for separation

*Old times:
Ethidium bromide*



+ *cost-effective*

 *mutagenic*




Traditional PCR

In traditional PCR the product is detected only after amplification is completed, using methods as gel electrophoresis





Gel electrophoresis separates non-target products from target based on size difference



Non-target amplification is commonly observed in gel-based PCR

PCR products visualization



Modern times:

Non-toxic dyes, such as GelRed, RedTaq, etc

WHAT IS PCR USED FOR?

Consumer genomics



Environmental microbiology



Medicine



Genetic research




Food and agriculture



Forensic science



Phylogenetics




Environmental biology:

Detection of microorganisms in the environment:

- PCR allows specific target species to be identified and quantified, even when very low numbers exist

e.g. searching for pathogens or indicator species such as coliforms in water supplies



Consumer genomics:


PCR has enabled personalized genome testing. An industry has sprung up offering consumers tailor-made products and services based on information in their genome. e.g. nutrigenomics



- **Nutrigenomics** is concerned with the impact of dietary components on the genome, the proteome (the sum total of all proteins), and the metabolome (the sum of all metabolites)
- Much of the nutrigenomic focus has been on single-nucleotide polymorphisms (SNPs): Dietary factors may differentially alter the effect of one or more SNPs to increase or decrease disease risk

SNPs: DNA sequence variations that account for 90% of all human genetic variation. SNPs that alter the function of “housekeeping genes” involved in the basic maintenance of the cell are assumed to alter the risk of developing a disease.

PCR applications



LIVE
EOT


Who's the father?

Mother

Child

F1

F2



NCIS

INSIDE THE LAB


Food and agriculture

Genetic technologies include a range of techniques that enable the modification of existing organisms for the purpose of improving foods and food production.

Although selective breeding has been around for centuries, with the advent of technology genetic manipulation can be more targeted and new crops produced more quickly.

e.g. transgenic cows: The transgene is present in every cell in the transgenic cow. However, it's only expressed in mammary tissue. This means that the transgene's protein will only be found in the cow's milk.

Increase the Yield of Milk




Genetic research

There are many fields of genetic research that use PCR as an essential tool. The majority of these technologies have multidisciplinary applications. These include the creation, detection and monitoring of genetically modified organisms (GMOs), genetic engineering, gene modification, transgenics, cloning, synthetic biology and directed evolution.

PCR in medicine:

- **Genetic disorders:** for the identification of many genes related to certain phenotypes, including genetic disorders. Genetic disorders are caused by mutations that range from simple changes in the base sequence of the DNA double helix through to changes in large DNA sequences and even whole chromosomes.
- The **Human Genome Project** for identifying and mapping all of the genes of the human genome
- **Genotyping** for sex determination of embryos as well as detecting chromosomal and genetic disorders in the fetus
- **Personalized medicine:** to tailor therapy with the best response and highest safety margin to ensure better patient care
- **Pharmaceutical production:** some medicinal products are produced by genetic modification (GM) of existing organisms. This includes GM microorganisms to mass produce pharmaceuticals such as insulin and some antibiotics.
- **Transgenics:** involves transfer of genes from one organism to another via genetic modification
- **Transplant testing:** technologies can use PCR to improve the accuracy of tissue compatibility, thereby reducing the likelihood of transplant rejection
- **Pathogenic microorganisms**, including some viruses, bacteria, parasites and fungi, cause infectious diseases and can be identified using PCR, aiding efficient diagnosis and treatment. e.g testing for Legionnaires' disease, tracing zoonoses and SARS-CoV-2

WHAT IS PCR USED FOR?




Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2) causes the COVID-19 disease.

PCR tests detect the genetic material of the virus (molecular diagnostics)



How to monitor the levels in gene expression?

Real-time QPCR measures the fluorescence at each cycle as the amplification progresses




Real-time PCR

Real-time PCR is a technology that allows amplification to be detected automatically at each cycle



Real-time Quantitative PCR


Real-time Quantitative PCR is a method to measure gene concentrations in samples



Required components for Real-Time PCR




Real-Time PCR Instrument:
Thermal cycler & optical detector



Real-Time PCR Chemistry:
Fluorescent signal automatically
emitted each cycle



Benefits of Real-Time PCR

- No laboratory procedures post-PCR
- More information than gels
- **Gene quantification with accuracy, precision and throughput**




PCR applications

WHAT IS PCR USED FOR?




Relative positions of amplicon targets on SARS-CoV-2 (previous Wuhan-CoV) genome. ORF: open reading frame; RdRp: RNA-dependent RNA polymerase; S: surface glycoprotein; N: nucleocapsid. Numbers below amplicon are genome locations corresponding to Wuhan-CoV genome, [MN908947].





Real-time PCR: SYBR Green Melting Curve

For SYBR Green assays where specificity is important, a melting curve should be performed after each amplification and the melting curve data evaluated


For example, a narrow symmetrical melting peak that overlays with a known pure target melting peak is consistent with single target amplification



Real-time PCR: SYBR Green Melting Curve (or Dissociation curve)




Real-time PCR: SYBR Green Melting Curve (or Dissociation curve)




Real-time PCR: SYBR Green Melting Curve

For SYBR Green assays intended to amplify one target, data from wells with multiple melt peaks should be considered compromised




Real-time PCR: 3 phases




Real-time PCR: Exponential phase

- Exponential (or geometric or logarithmic) phase is the first phase of real-time PCR where the efficiency is constant
- Quantitative data are derived from this phase




Real-time PCR: Linear phase

- Linear is the second phase where the efficiency continually declines, cycle to cycle



Real-time PCR: Plateau phase


- Plateau is the last phase where there is no significant amplification



Real-time PCR: reaction efficiency


The exponential phase has a unique behavior:
efficiency is not affected by initial target quantity

Reaction Efficiency is measured in the Exponential phase




Real-time PCR: threshold

Threshold is a horizontal line drawn through the exponential phase data




Threshold cycle (Ct or Cq) is the PCR cycle where the exponential data and threshold intersect




*Ct values are usually not whole numbers.

principles of real-time PCR: Ct & amplification plot



Ct: The first cycle at which the instrument can distinguish the amplification generated fluorescence as being above the ambient background


The greater the amount of initial DNA template in the sample, the earlier the Ct value for that sample



Real-time PCR: Ct & target quantity

Ct marks initial target quantity

A mathematical relationship exists between Ct and initial target quantity



$$\text{Quantity} \sim e^{-Ct}$$

e = geometric efficiency
e = 2 in this example.

Do I need a standard curve?




Figure 20
Amplification plots of standards in a four-fold dilution series over three orders of magnitude.

A standard curve enables you to identify the linear working range and the efficiency of your assay

- Perform a serial dilution series over a wide range of concentrations using a positive sample/control
- Use replicates for your standard curve
→ enables outlier detection and statistics

Properties of a good standard curve:


- high efficiency (90%-110%)
- good R^2 (>0.98)
- low replicate variability for individual standards



Since there is a *theoretical doubling* of the starting amount in each cycle what is the *efficiency*? $E \sim 2$

Real-time PCR: target concentration

- Quantification of gene targets can be defined as the concentration of the target in a biological sample
- Biological events in cells and tissues are based on changes in the concentration of molecules





Real-time PCR: absolute or relative quantification

- ❖ Real-time PCR experiments can be used to produce either relative or absolute data -> semi-quantitative or quantitative real-time PCR


Absolute quantification: molecular concentration of the target is calculated

Relative quantification: concentration of the target is calculated using a non-molecular or arbitrary unit

Application	Measurement Goal	Type
Gene Expression	Changes in mRNA concentration	Relative
Virus quantification	Molecular concentration of virus	Absolute
Transgenic copy #	Copies transgene/genome	Absolute



What Is Your Goal?



Comparative or Relative quantitation:

- ❖ comparing the levels of a target gene in test samples relative to a sample of reference

Quantitative as Absolute quantitation:

- a standard curve to quantitate the amount of target present in an Unknown sample
- or a series of Standard samples, containing dilutions of a known amount of target

The standard curve is then used to **derive the initial template quantity in Unknown wells based on their Ct values**

Real-time PCR: absolute or relative normalization

Absolute quantity (AQ) values are expressed as molecules of target in the numerator per normalization unit in the denominator




1300 molecules target/ng total RNA

Examples are:

- ng of total DNA
- copies of genomic DNA
- cell number

Relative quantity (RQ) values are expressed as numbers without a denominator



For example, a RQ of 2 means that the sample has twice the concentration of the target compared to a sample with a RQ of 1

Real-time PCR: normalizers

- A normalizer (or housekeeping or reference or control) is a gene used in real-time quantitative PCR experiment to normalize sample mass

Examples

18S rRNA
Acidic ribosomal protein
Beta-actin
Cyclophilin
Glyceraldehyde-3-phosphate dehydrogenase
Phosphoglycerokinase
β_2 -Microglobulin
β -Glucuronidase
Hypoxanthine ribosyl transferase
Transcription factor IID, TATA binding protein
Transferrin receptor

The expression levels of reference genes should remain constant between the cells of different tissues, the developmental stages, under different experimental conditions, etc.

- In order to compare samples within a group, normalization should be based on the same reference gene

Real-time PCR: endogenous and exogenous normalizers

Normalizers are categorized into two types:

- ✓ Endogenous
- ✓ Exogenous

- Endogenous normalizers are genes naturally present in each sample at a concentration relevant for normalization
- Endogenous normalizers are commonly used for samples derived from cells or tissues

Examples include 18S ribosomal protein to measure eukaryotic mass

- Exogenous normalizers are genes that are introduced to the sample at a concentration relevant for normalization
- Exogenous normalizers are typically used for liquid samples, such as blood, plasma, etc
- The exogenous normalizers are genes not naturally found in the samples being studied

Real-time PCR: normalizers' functions

Additional functions of the normalizer genes

- Correct for extraction efficiency of liquid samples
- Assist assay assessments and troubleshooting
- Amplification control
- Reverse transcription control
- Correct for pipette calibration error



Real-time PCR: normalization

- Normalization using quantity values is a division process
- Normalizing using Ct values is a subtraction process. It is used in the $\Delta\Delta C_t$ method

$$\frac{\text{Target Quantity}}{\text{Normalizer Quantity}}$$

$$\text{Target Ct} - \text{Normalizer Ct} = \Delta C_t$$

Real-time PCR: normalization

- Normalization using quantity values is a division process
- Normalizing using Ct values is a subtraction process. It is used in the $\Delta\Delta C_t$ method

$$\frac{\text{Target Quantity}}{\text{Normalizer Quantity}}$$

$$\text{Target Ct} - \text{Normalizer Ct} = \Delta C_t$$

Real-time PCR: normalization using the $\Delta\Delta C_t$ method

- $\Delta\Delta C_t$ is the recommended method to calculate real-time data
- $\Delta\Delta C_t$ is based on the equation: **Quantity = $2^{-\Delta\Delta C_t}$**

Sample	Detector	Task	Ct	ΔC_t	Avg ΔC_t	ΔC_t SD	$\Delta\Delta C_t$	RQ
Sample 2	B	Target	34.034378		1.07405205E1	1.7314075E-1	8.1801796E-1	5.6722069E-1
Sample 2	B	Target	34.711212		1.07405205E1	1.7314075E-1	8.1801796E-1	5.6722069E-1
Sample 2	H	Target	30.378319		6.514831543	1.9364803E-1	-1.3458252E-2	1.0093722343
Sample 2	H	Target	29.62249		6.514831543	1.9364803E-1	-1.3458252E-2	1.0093722343
Sample 2	V	Target	27.189629		3.0258388519	2.184637E-1	1.1824417E-1	9.2130822E-1
Sample 2	V	Target	26.394093		3.0258388519	2.184637E-1	1.1824417E-1	9.2130822E-1

Real-time PCR: normalization using the $\Delta\Delta C_t$ method

Normalisation according to an internal reference gene

“delta-delta Ct method” for comparing relative expression results between treatments in real-time PCR

ABI Prism Sequence detection System User Bulletin #2 (2001)

Relative quantification of gene expression

$$\Delta CP = CP_{\text{target gene}} - CP_{\text{reference gene}}$$

$$\text{expression ratio} = 2^{-[\Delta CP_{\text{treatment}} - \Delta CP_{\text{control}}]}$$

$$\text{expression ratio} = 2^{-\Delta\Delta CP}$$

Livak KJ, Schmittgen TD. (2001)

Analysis of Relative Gene Expression Data Using Real-Time Quantitative PCR and the $2^{-\Delta\Delta C(T)}$ method.

Methods, 2001 25(4): 402-408.



Important for the analysis:

Include a negative control & a positive control

PCR controls

Positive controls

- Necessary for initial validation
- help to generate performance data: linear range, sensitivity

Examples:

PCR products, linearized plasmids, positive samples

Negative controls

- Necessary for initial validation and troubleshooting
- address specificity

Examples:

No Template Control (NTC)
→ Primer dimers, probe degradation

No RT Control (NoRT)
→ Genomic amplification

Controls should identify which data is good or bad!

How can the products be spotted?


SYBR-Green chemistry:

- cost-effective
- easy to use

Unbound SYBR® Green




Fluorescent Bound SYBR® Green



Probe-Based chemistries (e.g. Taqman)

- **higher detection specificity**
- multiple probes labeled with different reporter dyes allow detection of more than one target in a single reaction (**multiplex QPCR**)




The fluorescence intensity increases proportionally with each amplification cycle in response to the increased amplicon concentration



Important for the analysis:

Sequences of primers or probes should be removed!

qRT-PCR reactions: One-step vs Two-step: limiting sample amount?



Design an experiment!



To considerate:

1. is this an absolute or relative quantification?
2. what is the target gene and what is the reference gene?
3. how many samples
4. how many replicates? Monoplicates / duplicates / triplicates, etc for each sample?
5. positive / negative controls? If yes, how many?
6. how will you choose the appropriate cycling conditions? What will they be?

INVESTIGATOR:

RUN No:

DATE:

	1	2	3	4	5	6	7	8	9	10	11	12
A	○	○	○	○	○	○	○	○	○	○	○	○
B	○	○	○	○	○	○	○	○	○	○	○	○
C	○	○	○	○	○	○	○	○	○	○	○	○
D	○	○	○	○	○	○	○	○	○	○	○	○
E	○	○	○	○	○	○	○	○	○	○	○	○
F	○	○	○	○	○	○	○	○	○	○	○	○
G	○	○	○	○	○	○	○	○	○	○	○	○
H	○	○	○	○	○	○	○	○	○	○	○	○

Copyright © 2004 Inno Assenboort

Samples:

3 patients with breast cancer: P1, P2, P3

Reagents:

Reaction volumes per well:	
cdDNA	3uL
SYBR Green Mix	17uL
Forward primer (10mM)	0.24uL
Reverse primer (10mM)	0.24uL
ddH ₂ O	up to 20mL

mRNA targets:

- 1) BRAF
- 2) GADPH

Results:

Design an experiment!



To considerate:

is this an absolute or relative quantification?

relative quantification

what is the target gene and what is the reference gene?

target gene: BRAF, reference gene: GADPH

how many samples & how many replicates? Monoplicates / duplicates / triplicates, etc. for each sample?

it depends!

positive / negative controls? If yes, how many?

negative controls: NTC, NoRT, positive controls: any gene that you know it is stably expressed

how will you choose the appropriate cycling conditions? What will they be?

that depends on SYBR Green qPCR Kit Master Mix

INVESTIGATOR:

RUN No:

DATE:

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

Samples:

3 patients with breast cancer: P1, P2, P3

Reagents:

Reaction volumes per well:

cdNA	3uL
SYBR Green Mix	17uL
Forward primer (10mM)	0.24uL
Reverse primer (10mM)	0.24uL
ddH ₂ O	up to 20mL

mRNA targets:

- 1) BRAF
- 2) GADPH


Results:

microfluidics

The principle


When a particle passes through a laser beam, it scatters light according to its physical properties.

e.g. A larger particle will cause higher deviation of the light. If this particle is also fluorescent, after being excited with light of the appropriate wavelength, it will emit light at another (usually higher) wavelength.



For cells

- Cells should be in single-cell suspension
- Cells suspended in a liquid (fluidics system) pass through a laser beam causing light scatter that can be sensed by optical-to-electronic coupling system detectors
- The electronics system process the signal from the detectors giving data that represent the features (size and granularity) of the cells



microfluidic PCR

Microfluidic devices offer small thermal mass, low thermal inertia, and rapid heat transfer:

- The **small volumes reduce sample and reagent consumption**, leading to inexpensive operation of the system
- More **efficient mass transport and heat conduction** in microfluidic devices: while in macroscale systems diffusion and heat conduction are slow for efficient transport, the smaller scales of microfluidic devices allow for utilization of these effects

Microfluidics have been used for PCR, qPCR, RT-PCR

Lab-on-a-Chip


microfluidics benefits: large number of samples, speed, parallelization, sensitivity, multiplexity, reduction of reagents cost

microfluidics drawbacks: PCR inhibition and carryover contamination

microfluidic PCR examples

Various designs and approaches of microfluidic PCR devices

Radial PCR
water-in-oil emulsion



85bp DNA are amplified in 17min within 34 cycles

denaturation

annealing extension

microfluidic PCR: Droplet Digital PCR (ddPCR)

Endless Frontier of 'Divide and Conquer'

- Droplet Digital PCR (ddPCR) is a method for performing PCR that is based on **water-oil emulsion** droplet technology
- ddPCR technology uses a combination of microfluidics and proprietary surfactant chemistries to divide PCR samples into water-in-oil droplets. **A sample is fractionated into thousands of nanoliter-sized droplets, and PCR amplification of the template molecules occurs in each individual droplet**
- This partitioning enables the measurement of thousands of **independent amplification events within a single sample**
- The droplets serve essentially the same function as individual test tubes or wells in a plate in which the PCR reaction takes place, albeit in a much smaller format



microfluidic PCR: Droplet Digital PCR (ddPCR)

- Following PCR, **each droplet is analyzed or read to determine the fraction of PCR-positive droplets in the original sample**. These data are then analyzed using Poisson statistics to determine the **target DNA template concentration** in the original sample
- ddPCR technology uses reagents and workflows similar to those used for most standard TaqMan probe-based assays. This technique has a **smaller sample requirement** than other commercially available digital PCR systems, reducing cost and preserving precious samples



One measurement



Nanodroplet PCR reactions are independent, single amplification events




Many thousands of discrete measurements

microfluidic PCR: Droplet Digital PCR (ddPCR)

DNA undergoes limiting dilution

- partition volumes: >5 pL
- the number of partitions: >100.000 for a single experiment



- 1 Isolate nucleic acid and quantitate
- 2 Mix samples with real-time PCR reagents and prepare dilutions
- 3 Split reaction mixture into many partitions
- 4 Amplify and collect data
- 5 Analyze data to identify copies/reaction



microfluidic PCR: Droplet Digital PCR (ddPCR)

ddPCR should be performed when

- ☺ rare events must be quantified
- ☺ individual mutations must be detected

ADVANTAGES of ddPCR:

- No need to rely on references or standards
- Absolute quantification
- Increased precision
- Removal of PCR bias
- Capability to analyze complex mixtures
- Linear detection of small-fold changes

ddPCR APPLICATIONS:

- ✓ viral detection
- ✓ diagnostic monitoring
- ✓ fetal screening
- ✓ *research applications* such as phage-host interactions and intracellular profiling

ddPCR data output




Fig. 1. Examples of data output from a droplet ddPCR instrument (Bio-Rad QX100). Clear distinction between **positive and negative partitions (A and B)**; **no-template control reactions (C and D)**.

organ-on-a-chip


An **organ-on-a-chip** is a **microfluidic cell culture device** created with microchip manufacturing methods that contains continuously perfused chambers inhabited by living cells arranged to simulate tissue- and organ-level physiology.

By recapitulating the multicellular architectures, tissue-tissue interfaces, physicochemical microenvironments and vascular perfusion of the body, these devices produce levels of tissue and organ functionality not possible with conventional 2D or 3D culture systems.

They also enable high-resolution, real-time imaging and *in vitro* analysis of biochemical, genetic and metabolic activities of living cells in a functional tissue and organ context.






Applications

- Study tissue development, organ physiology and disease etiology
- Study molecular mechanisms of action, toxicity testing and biomarker identification in the area of drug discovery and development



lung-on-a-chip

A microfluidic chip is suitable for growing lung-on-a-chip cell cultures, with an upper channel used for air flow, a lower channel used for liquid flow, and side chambers to apply different pressures



(d) A ‘breathing’ lung-on-a-chip that recapitulates the alveolar-capillary interface. Human alveolar epithelial cells are cultured on top of a flexible, porous, ECM-coated membrane and human capillary endothelial cells on the bottom. Air is passed through the upper channel to create an air-liquid interface with the alveolar epithelium, and culture medium is flowed through the vascular channel, with or without human immune cells (top). Breathing motions are mimicked by applying cyclic suction to full-height side chambers that rhythmically distort and relax the flexible side walls and attached porous membrane. The fluorescence confocal 3D reconstruction at the bottom and the cross-sectional view shown at higher magnification in the inset show the tissue-tissue interface formed between the alveolar epithelium (green) and the endothelium (red).

PCR algorithms, packages...

for housekeeping genes:

geNorm, NormFinder, BestFinder, BestKeeper, Geninvestigator, RefGenes, REST

R packages:

SLqPCR, EasyqpcR, ddCt, NormqPCR, ReadqPCR, FPK-PCR, HTqPCR

For extra info

[Science](#). 1991 Jun 21;252(5013):1643-51.

Recent advances in the polymerase chain reaction.

[Erlich HA](#)¹, [Gelfand D](#), [Sninsky JJ](#).

[Lab Chip](#). 2016 Oct 5;16(20):3866-3884.

Polymerase chain reaction in microfluidic devices.

[Ahrberg CD](#)¹, [Manz A](#)², [Chung BG](#)¹

[Methods Mol Biol](#). 2014;1169:133-42. doi:

10.1007/978-1-4939-0882-0_13.

Expression profiling by real-time quantitative polymerase chain reaction (RT-qPCR).

[Lech M](#)¹, [Anders HJ](#).

[Expert Rev Mol Diagn](#). 2005 Mar;5(2):209-19.

Basic principles of real-time quantitative PCR.

[Arya M](#)¹, [Shergill IS](#), [Williamson M](#), [Gommersall L](#),
[Arya N](#), [Patel HR](#).

[Nat Protoc](#). 2006;1(3):1559-82.

Quantification of mRNA using real-time RT-PCR.

[Nolan T](#)¹, [Hands RE](#), [Bustin SA](#).

[Anal Biochem](#). 2005 Sep 1;344(1):141-3.

The implications of using an inappropriate reference gene for real-time reverse transcription PCR data normalization.

[Dheda K](#)¹, [Huggett JF](#), [Chang JS](#), [Kim LU](#), [Bustin SA](#), [Johnson MA](#),
[Rook GA](#), [Zumla A](#).

[Methods](#). 2001 Dec;25(4):402-8.

Analysis of relative gene expression data using real-time quantitative PCR and the 2(-Delta Delta C(T)) Method.

[Livak KJ](#)¹, [Schmittgen TD](#).

[Nucleic Acids Res](#). 2001 May 1;29(9):e45.

A new mathematical model for relative quantification in real-time RT-PCR.

[Pfaffl MW](#)¹.

For extra info

[J Vis Exp.](#) 2012 May 22;(63):e3998. doi: 10.3791/3998.

Polymerase chain reaction: basic protocol plus troubleshooting and optimization strategies.

[Lorenz TC](#)¹.

[Sci Rep.](#) 2017 Aug 29;7(1):9617. doi: 10.1038/s41598-017-09183-4.

Poisson Plus Quantification for Digital PCR Systems.

[Majumdar N](#)¹, [Banerjee S](#)², [Pallas M](#)², [Wessel T](#)², [Hegerich P](#)².

[Clin Chem.](#) 2013 Jun;59(6):892-902. doi:

10.1373/clinchem.2013.206375. Epub 2013 Apr 9.

The digital MIQE guidelines: Minimum Information for Publication of Quantitative Digital PCR Experiments.

[Huggett JF](#)¹, [Foy CA](#), [Benes V](#), [Emslie K](#), [Garson JA](#), [Haynes R](#), [Hellemans J](#), [Kubista M](#), [Mueller RD](#), [Nolan T](#), [Pfaffl MW](#), [Shiple GL](#), [Vandesompele J](#), [Wittwer CT](#), [Bustin SA](#).

[Clin Chem.](#) 2015 Jan;61(1):79-88. doi:

10.1373/clinchem.2014.221366. Epub 2014 Oct 22.

Considerations for digital PCR as an accurate molecular diagnostic tool.

[Huggett JF](#)¹, [Cowen S](#)², [Foy CA](#)².

[Anal Chem.](#) 2011 Nov 15;83(22):8604-10. doi:

10.1021/ac202028g. Epub 2011 Oct 28.

High-throughput droplet digital PCR system for absolute quantitation of DNA copy number.

[Hindson BJ](#)¹, [Ness KD](#), [Masquelier DA](#), [Belgrader P](#), [Heredia NJ](#), [Makarewicz AJ](#), [Bright IJ](#), [Lucero MY](#), [Hiddessen AL](#), [Legler TC](#), [Kitano TK](#), [Hodel MR](#), [Petersen JF](#), [Wyatt PW](#), [Steenblock ER](#), [Shah PH](#), [Bousse LJ](#), [Troup CB](#), [Mellen JC](#), [Wittmann DK](#), [Erndt NG](#), [Cauley TH](#), [Koehler RT](#), [So AP](#), [Dube S](#), [Rose KA](#), [Montesclaros L](#), [Wang S](#), [Stumbo DP](#), [Hodges SP](#), [Romine S](#), [Milanovich FP](#), [White HE](#), [Regan JF](#), [Karlin-Neumann GA](#), [Hindson CM](#), [Saxonov S](#), [Colston BW](#).

[Anal Chim Acta.](#) 2016 Mar 31;914:7-16. doi:

10.1016/j.aca.2016.02.006. Epub 2016 Feb 13.

A review on continuous-flow microfluidic PCR in droplets: Advances, challenges and future.

[Zhang Y](#)¹, [Jiang HR](#)².

[Lab Chip.](#) 2016 Oct 5;16(20):3866-3884.


Polymerase chain reaction in microfluidic devices.

[Ahrberg CD](#)¹, [Manz A](#)², [Chung BG](#)¹.

Polymerase Chain Reaction (PCR)

key-notes

PCR steps: 1. denaturation, 2. annealing, 3. extension

PCR 

- Absolute with a standard curve
- Relative with reference genes

endpoint PCR
digital PCR
microfluidic PCR

VS.

Real-time PCR

Polymerase Chain Reaction (PCR)



questions?

ТНАОК



hands on!

Learn



unity


A/a

Design

HARE

* Focus on ... ty!


At the bench side..



CGCCCAATACGACCAAATCC

Real-time PCR: Standard curve for dynamic range testing

- Dynamic range is the minimum and maximum sample amount that may be run in a quantitative assay and produce accurate and precise data
- Before starting a new quantitative experiment, dynamic range should be tested
- Dynamic range testing involves standard curves



Real-Time PCR applications -> SNP genotyping

- **Single Nucleotide Polymorphisms (SNPs)** are single base changes in genomic DNA
- **SNPs may cause disease or regulate metabolic traits**
- **SNPs have consisted prominent biomarkers**, e.g. in patients with localized renal cell cancer shorter overall survival (VEGFR2: rs10013228; VEGFR3: rs6877011) ($p < 0.05$)

SNPs in diploid loci have 3 genotypes:

Homozygous

A

A

Heterozygous

A

G

Homozygous

G

G

Real-Time PCR applications -> Quantitative applications

Quantitative applications:

- Gene expression
- Transgenic copy number
- Cancer
- Forencics
- GMO testing
- Quality control
- Pathogen quantification
- Chimerism
- Next Gen library quantification



Real-time PCR -> Allelic discrimination

- Allelic discrimination is a real-time PCR chemistry for SNP genotyping
- Allelic discrimination uses two probes: one for each allele
- Analysis is performed on post-PCR read results

